VitriFreeze EST™ / VitriThaw EST™

INTENDED USE

VitriFreeze EST™ and VitriThaw ES™ are a set of media for vitrification and thawing of human embryos (zygote till expanded blastocyst stage). For professional use only.

MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for a maximum of 6 vitrification/warming procedures, depending on the protocol used.

VitriFreeze ES™ kit (product code VF_KIT1_ES): contains one bottle of the following media:

- Product code VPIE005: 5 ml VitriFreeze ES Pre-incubation medium ("PI")
- Product code VTE001: 1 ml VitriFreeze ES Freezing medium 1 ("F1")
- Product code VFE2001: 1 ml VitriFreeze ES Freezing medium 2 ("F2")
- Product code VFE3001: 1 ml of Vitrifreeze ES Freezing medium 3 ("F3")

VitriThaw ES™ kit (product code VT_KIT1_ES): contains one bottle of the following media:

- Product code VTE1005: 5 ml of VitriThaw ES Thawing medium 1 ("T1")
- Product code VTE2001: 3.2 ml of VitriThaw ES Thawing medium 2 ("T2")
- Product code VTE3001: 1 ml of VitriThaw ES Thawing medium 3 ("T3")
- Product code VTE4001: 1 ml of VitriThaw ES Thawing medium 4 ("T4")
- Product code VTE5001: 1 ml of VitriThaw ES Thawing medium 5 ("T5")

The media should be used in the order displayed above (the bottles may be in a different order in the kit).

MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device, preferably closed system (VitriFreeze ES™/VitriThaw ES™ were clinically validated with use of the following closed devices: High Security Vitrification straw (HSV, CBS) or Vitrisafe (Nextclinics). Importantly, each laboratory should validate their own procedures when using VitriFreeze ES™/VitriThaw ES™ in combination with other vitrification devices)

COMPOSITION

VitriFreeze ES™ and VitriThaw ES™ are ready to use. The media contain PBS and human serum albumin (11-20 g/liter depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening for individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens. Aseptic techniques should be used to avoid possible contamination.

Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

PRODUCT SPECIFICATIONS

- Chemical composition
 - pH: 7.20 - 7.40
 - Osmolality (mOsm/kg):
 - VitriFreeze ES Pre-incubation medium: 270-295 (release criteria: 270-290)
 - VitriThaw ES Thawing medium 2: 805-865 (release criteria: 805-850)
 - VitriThaw ES Thawing medium 3: 535-565
 - VitriThaw ES Thawing medium 4: 405-435
 - Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: no growth of endotoxins test by Limulus Lysate (LAL) methodology (USP <85>): <0.25 EU/ml
 - Mouse Embryo Assay (one-cell mouse embryo development until blastocyst stage after 96 hours of exposure to the cooling and warming media) ≥ 80%
 - Use of Ph Eur or USP grade products if applicable
 - The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for VitriFreeze ES™ and VitriThaw ES™ are ready to use. The media contain PBS and human serum albumin (11-20 g/liter depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

VitriFreeze ES™ media have increasing concentrations of DMSO and Ethylene glycol (EG). VitriFreeze ES™ Freezing 3 also contains sucrose and ficol. VitriThaw ES™ media have decreasing concentrations of sucrose.

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- Stocker entre 2-8°C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm the media to room temperature (20-25°C).

• Open the necessary number of vitrification devices, taking into account that 1 device can hold 1-2 embryos, (for maximum volume load: check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.

• Prepare the following dishes for vitrification or warming of embryos as described in figure 1.

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Notes:

- Do not use the same droplets/ media set-up for different patients.

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