VitriFreeze EST™ / VitriThaw EST™

INTENDED USE

VitriFreeze EST™ / VitriThaw EST™

MATERIAL INCLUDED WITH THE KIT

- One kit provides sufficient medium for a maximum of 6 vitrification/warming procedures, depending on the protocol used.
- VitriFreeze EST™ kit (product code VF_KIT1_ES): contains one bottle of the following media:
 - Product code VPIE005: 5 ml VitriFreeze ES Pre-incubation medium ("PI")
 - Product code VTE2001: 3.2 ml of VitriThaw ES Freezing medium 1 ("F1")
 - Product code VTE3001: 1 ml VitriFreeze ES Freezing medium 2 ("F2")
 - Product code VFE3001: 1 ml of VitriFreeze ES Freezing medium 3 ("F3")
- VitriThaw EST™ kit (product code VT_KIT1_ES): contains one bottle of the following media:
 - Product code VTE1005: 5 ml VitriThaw ES Thawing medium 1 ("T1")
 - Product code VTE2001: 3.2 ml of VitriThaw ES Thawing medium 2 ("T2")
 - Product code VTE3001: 1 ml of VitriThaw ES Thawing medium 3 ("T3")
 - Product code VTE4001: 1 ml of VitriThaw ES Thawing medium 4 ("T4")
 - Product code VTE5001: 1 ml of VitriThaw ES Thawing medium 5 ("T5")

The media should be used in the order displayed above (the bottles may be in a different order in the kit).

MATERIAL NOT INCLUDED BUT NEEDED

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. Known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard laboratory safety and infection control practices, including the use of medical products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

COMPOSITION

VitriFreeze EST™ and VitriThaw EST™ are ready to use. The media contain PBS and human serum albumin (11.20 g/liter) depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

VitriFreeze EST™ media have increasing concentrations of DMSO and Ethylene glycol (EG). VitriFreeze ES Freezing 3 also contains sucrose and ficol. VitriThaw EST™ media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
pH: 7.20 - 7.40
- Osmolality (mOsm/kg):
• VitriFreeze ES Pre-incubation medium: 270-295 (release criteria: 270-290)
- VitriThaw ES Thawing medium 2: 805-865 (release criteria: 805-850)
- VitriThaw ES Thawing medium 3: 535-565
- VitriThaw ES Thawing medium 4: 405-435
- Sterility test by the current Ph. Eur. 2.6.1./ USP <71>: no growth
- Endotoxins test by Limulus Lysate (LAL) methodology (USP <85>): <0.25 EU/ml
- Mouse Embryo Assay (one-cell mouse embryo dechorionated blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

STERILE A

Document ID: FP09 I46 02 R01 E.2

Update: 21/03/2023

MATERIAL INCLUDED / MATERIEL INCLUS / IM LIEFERUMFANG ENTHALTENE MATERIALIEN / MATERIAL INCLUIDO / MATERIALI INCLUSI / MATERIAL INCLUIDO / YUKO ΠΟΥ ΠΕΡΙΛΑΜΒΑΝΕΤΑΙ / MATERIAAL AANWEZIG

Catalogue number	Référence catalogue	Número de catálogo	Numero di catalogo	Aphelyt-katalognummer	Catalogusnummer
VF_KIT1_ES	VitriFreeze ES - kit				
VT_KIT1_ES	VitriThaw ES - kit				

CUSTOMER-TECHNICAL SUPPORT

SUPPORT CLIENTS-SUPPORT TECHNIQUE
KUNDENDIENST-TECHNISCHE SUPPORT
ATENCIÓN AL CLIENTE-ASISTENCIA TÉCNICA
ASSISTENZA CLIENTI-SUPPORTO TECNICO
APOYO AL CLIENTE-TECNICO
ΥΠΟΣΤΗΡΙΞΗ ΠΕΛΑΤΩΝ-ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ
KLANTEN-TECHNISCHE ONDERSTEUNING

- Notes:**
- Do not use the same droplets/ media set-up for different patients.
 - Set-up with 4-well dish (for cases with multiple embryos): use 5 vitrification cycles of the same patient can be performed in the same media.
 - Set-up with droplets (for cases with few embryos): the number of embryos that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.

Vitrifreeze protocol

- Expose the embryos sequentially to the following media as described in [figure 2](#).

- (*) Recommended pipetting technique for step using F3 ([figure 3](#)):

- Place the embryo on the bottom of the 1st drop of F3 (step A); Make sure the embryo is in the drop. Preferably, do not extrude the complete volume of the F2. The embryo will float.
- Empty the pipette outside the drop (step B1).
- Rinse the pipette with medium from drop of F2 (step B2).
- Aspirate a small volume of drop 2, aspirate the embryo from drop 1, transfer to the bottom of drop 2 (step C).
- Make sure the embryo is in the drop. Flush in and out 2 times.
- Place the embryo in drop 3 and flush in and out 3 times until the shadow-like shape disappears (step D).
- Now load the embryo on the vitrification device.

- Thawing protocol

- Expose the embryos sequentially to the following media as defined in [figure 4](#).

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. Known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

- Standard laboratory safety and infection control practices, including the use of medical products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

COMPOSITION

VitriFreeze EST™ and VitriThaw EST™ are ready to use. The media contain PBS and human serum albumin (11.20 g/liter) depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

VitriFreeze EST™ media have increasing concentrations of DMSO and Ethylene glycol (EG). VitriFreeze ES Freezing 3 also contains sucrose and ficol. VitriThaw EST™ media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
pH: 7.20 - 7.40
- Osmolality (mOsm/kg):
• VitriFreeze ES Pre-incubation medium: 270-295 (release criteria: 270-290)
- VitriThaw ES Thawing medium 2: 805-865 (release criteria: 805-850)
- VitriThaw ES Thawing medium 3: 535-565
- VitriThaw ES Thawing medium 4: 405-435
- Sterility test by the current Ph. Eur. 2.6.1./ USP <71>: no growth
- Endotoxins test by Limulus Lysate (LAL) methodology (USP <85>): <0.25 EU/ml
- Mouse Embryo Assay (one-cell mouse embryo dechorionated blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

- The SSCP for VitriFreeze EST™ and VitriThaw EST™ describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by using the QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

WARNINGS BEFORE USE

MATERIAL INCLUDED DANS LE KIT

- Do not use the product if:
• it becomes cloudy, or shows any evidence of microbial contamination.
• seal of the container is opened or defect when the product is delivered.
• expiry date has been exceeded.
• Do not freeze before use.
• Do not re-sterilize after opening.
• Freezing medium 2 may contain small salt precipitates which do not have an impact on product performance/safety.

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STORAGE / DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from sunlight.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm the media to room temperature (20-25 °C).
- Open the necessary number of vitrification devices, taking into account that 1 device can hold 1-2 embryos, (for maximum volume load: check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of embryos as described in [figure 1](#).

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