

FertiCult™ High Viscosity Oil

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
LT	puslapis	11
TR	Sayfa	12
HU	oldal	13

FertiCult™ High Viscosity Oil

STERILE A

Document ID: FP09 I20V R01 B.2
Update: 08/05/2023

MATERIAL INCLUDED

Catalogue number		
HVOIL050P	1 x 50ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL100P	1 x 100ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL500P	1 x 500ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL050	1 x 50ml	FertiCult™ High Viscosity Oil (glass bottle)
HVOIL100	1 x 100ml	FertiCult™ High Viscosity Oil (glass bottle)
HVOIL500	1 x 500ml	FertiCult™ High Viscosity Oil (glass bottle)

CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

INDICATIONS FOR USE

FertiCult High Viscosity Oil is pre-washed with water, ready to use, high viscosity mineral oil (d=0.86-0.88g/ml) for use in In Vitro Fertilization (IVF), Intracytoplasmic Sperm Injection (ICSI) and Assisted reproductive technologies (ART).

Oil overlay is highly recommended during oocyte culture, denudation of oocytes, fertilization of oocytes (either by IVF or ICSI), embryo culture and manipulation. The oil permits small volumes of medium to be used for fertilization and embryo culture / handling. The primary function of an oil overlay is to prevent the culture medium from evaporation, which would result in increased osmotic pressure, and to avoid fluctuations in the pH and temperature. This fact is critical since variations in the pH of the medium may affect embryo metabolism, whereas osmotic stress may cause DNA and protein alterations that affect cell function.

For professional use only.

COMPOSITION

Mineral oil (Paraffin oil) (EP Monograph 0239, CAS code: 8042-47-5)

PRODUCT SPECIFICATIONS

- Chemical composition
- Density: 0.86-0.88g/ml
- Viscosity: 95-140cP at 30 °C
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85): < 0.1EU/ml
- Peroxide Value (POV): < 0.1mEq/kg (release criteria: <0.06 mEq/kg)
- Sterility test according to current Ph. Eur. 2.6.1/USP<71: no growth
- Mouse embryo assay: ≥ 80% blastocysts after 120 hours culture in standard mouse IVF medium covered with test item
- Use of pharmaceutical grade mineral (paraffin) oil
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use

STORAGE/DISPOSAL INSTRUCTIONS

- Store between 15-25 °C.
- Keep away from (sun)light and avoid exposure to sources of irradiation
- The products can be used up to 1 month after opening, when sterile conditions are maintained and the products are stored at room temperature.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

- FertiCult High Viscosity Oil is pre-washed. No further preparations are necessary.
- Compatible with gamete/embryo culture or handling media classified as class Ib or III Medical Device (for example but not limited to: FertiCult™ IVF medium, GAIN™ medium and other products of product range of ART media manufactured by FertiPro NV)
- We advise to prepare dishes/flasks with oil the day before use. Incubate them at 37 °C in a CO₂ incubator so that the oil warms up and becomes saturated with gas. If flasks are used, be sure not to tighten the cap of the flask because this will prevent gas exchange.
- Importantly: each laboratory should consult its own optimized and validated laboratory procedures.

PRECAUTIONS

- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens
- Aseptic techniques should be used to avoid possible contamination.
- The oil may appear cloudy if removed from the incubator and allowed to stand at room temperature for extended periods of time. Such cloudiness usually disappears if the oil is re-equilibrated by replacing it in the incubator for several hours, if it does not disappear discard the product.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

CLINICAL BENEFIT

The clinical benefit of FertiCult High Viscosity Oil is related to the fact that the oil is used during ART-procedures of patients and couples with infertility problems. It is expected that the outcomes of ART procedures in which the oil is used are comparable with the ART outcomes (i.e. pregnancy rate and delivery rate) as published in the annual peer-reviewed report of the European IVF-monitoring Consortium (EIM) for the European Society of Human Reproduction and Embryology (ESHRE).

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1

 Catalogue number

 Batch code

 Use-by date

 Manufacturer

 Keep away from sunlight

 Medical device

 Consult instructions for use

 Cautions

 Date of manufacture

 Sterilized using aseptic processing techniques

 Temperature limit

 Single sterile barrier system

 Single sterile barrier system with protective packaging outside.

 Do not use if package is damaged and consult instructions for use.


Symbol as defined in EU Commission MDR 2017/745

 CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

 US Federal law restricts this device to sale by or on the order of a physician.

FertiCult™ High Viscosity Oil

STERILE A

Document ID: FP09 I20V R01 B.2
Update: 08/05/2023

MATÉRIEL INCLUS

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HVOIL050P	1 x 50ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL100P	1 x 100ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL500P	1 x 500ml	FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL050	1 x 50ml	FertiCult™ High Viscosity Oil (glass bottle)
HVOIL100	1 x 100ml	FertiCult™ High Viscosity Oil (glass bottle)
HVOIL500	1 x 500ml	FertiCult™ High Viscosity Oil (glass bottle)

SUPPORT CLIENTS-SUPPORT TECHNIQUE

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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