FertiCult™ Aspiration medium

FertiCult™ Aspiration medium

FertiCult™ Aspiration medium (EN)

Cell culture medium
for human oocyte pick-up

STERILE A

Sterilised by sterile filtration.
Document reference: FP09 I19 R01 C.5
Update: 06.07.2020

GENERAL INFORMATION AND INTENDED USE

FertiCult™ Aspiration medium is a ready to use formulation for oocyte pick-up and/or follicular flushing.
FertiCult™ Aspiration medium contains HEPES; **no CO₂ incubation** is required.

As with all IVF media, FertiCult™ Aspiration medium has to be **pre-incubated in the incubator for 12 hours before use** (with lid closed).

For professional use only.

COMPOSITION

FertiCult™ Aspiration medium is a ready-to-use HEPES-buffered medium which contains bicarbonate, physiologic salts, glucose, lactate and 2.5 IU/ml heparin.

MATERIAL INCLUDED IN THE KIT

Product Code: ASPIO20
» 5x 20ml FertiCult™ Aspiration medium

Product Code: ASPH100
» 5x 100ml FertiCult™ Aspiration medium

MATERIAL NOT INCLUDED IN THE KIT

- » Incubator at 37°C (no CO₂)
- » Petri dishes
- » Microscope
- » Test tubes
- » LAF bench (ISO 5 environment)
- » Rubber free syringe (e.g. Braun Injekt)

PRODUCT SPECIFICATIONS

- » Chemical composition
- » pH: between 7.30-7.60 (Release criteria)
- » Osmolality: 270-290 mOsm/kg
- » Sterility: sterile (SAL 10⁻³)
- » Endotoxins: < 0.25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h):
≥ 80% after 60min of exposure (zygote stage)
- » Use of Ph Eur or USP grade products if applicable
- » Certificate of analysis and MSDS are available upon request

PRE-USE CHECKS

- » Do not use the product if it becomes cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

Symbols as defined in ISO 15223-1:2016
Symboles tels que définis dans ISO 15223-1:2016
Symbole nach der Definition von ISO 15223-1:2016
Simboli di cui alla norma ISO 15223-1:2016
Σμβολοs definidos en la norma ISO 15223-1:2016
Simbolos conforme definido na norma ISO 15223-1:2016
Symboelen volgens ISO 15223-1:2016
Σμβολα όπως ορίζονται στο ISO 15223-1:2016

STORAGE INSTRUCTIONS

- » Suitable for transport or short term storage at elevated temperatures (up to 5 days at 37°C)
- » Store between 2-25°C, once opened store between 2-8°C
- » Do not freeze before use
- » Keep away from (sun)light
- » The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C
- » Do not use after expiry date
- » Content cannot be re-sterilized after opening

WARNINGS AND PRECAUTIONS

All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination. FertiCult™ Aspiration medium contains heparin which is derived from porcine intestinal mucosa.

METHOD

OOCYTE COLLECTION

At oocyte collection, keep FertiCult™ Aspiration medium at 37°C or as close to this as possible.

Take about 20-50ml at a time, using a rubber free syringe (e.g. Braun Injekt). Assemble needle and tubing set according to manufacturer's instructions.

Attach syringe containing FertiCult™ Aspiration medium to the appropriate Luer, allowing connection to the lumen of the needle.

The person holding the syringe should try to apply warmth to it by holding it in the palm of the hand covering as much of the syringes surface as possible.

During aspiration of follicles, if an ovum is not extracted after a follicle has been emptied, flushing is indicated. The collapsed follicle may be refilled with FertiCult™ Aspiration medium, after which it may be aspirated once again. This may be repeated again and again, after refilling the syringe with fresh warm aspiration medium each time, until an ovum is recovered or the surgeon decides to abandon the follicle (after three attempts at flushing it may be sensible to move on, but the ultimate decision remains with the surgeon).

After oocyte collection, transfer oocytes to an appropriate culture medium (for example GAIN Medium, FertiPro N.V.).

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