EN FR DE IT

FertiCult™ IVF medium

FertiCult™ IVF medium

FertiCult™ IVF medium

STERILE A

Document ID: FP09 I07 R01 E.2
Update: 03.12.2019

USED ABBREVIATIONS

ICSI IntraCytoplasmatic Sperm Injection
IVF In Vitro Fertilization

GENERAL INFORMATION / INTENDED USE

FertiCult™ IVF medium is a bicarbonate buffered cell culture medium for washing and holding of human oocytes, performing fertilization by IVF or ICSI (until 2PN). The medium is complete and needs no further additives.

As with all IVF media, FertiCult™ IVF medium has to be pre-incubated in a CO₂ incubator with 6% CO₂ for at least 4 hours, but ideally, dishes should be prepared the day before use and incubated overnight (with lid opened) to obtain an optimal pH for oocyte washing/holding.

For professional use only.

COMPOSITION OF THE MEDIUM

FertiCult™ IVF medium is a ready-to-use bicarbonate-buffered balanced salt solution, supplemented with glucose (5.55 mM), lactate, pyruvate, and human serum albumin (4g/liter). The product is also available with phenol red and/or gentamicin (10mg/liter).

MATERIALS INCLUDED

Product code: FECU020

- » 5x 20ml FertiCult™ IVF medium

Product code: FECU050

- » 5x 50ml FertiCult™ IVF medium

Product code: FECU100

- » 3x 100ml FertiCult™ IVF medium

Product code: FECU020PHR

- » 5x 20ml FertiCult™ IVF medium with phenol red

Product code: FECU050PHR

- » 5x 50ml FertiCult™ IVF medium with phenol red

Product code: FECU100PHR

- » 3x 100ml FertiCult™ IVF medium with phenol red

Product code: FECU020PHR_G

- » 5x 20ml FertiCult™ IVF medium with phenol red and gentamicin

Product code: FECU050PHR_G

- » 5x 50ml FertiCult™ IVF medium with phenol red and gentamicin

Product code: FECU100PHR_G

- » 3x 100ml FertiCult™ IVF medium with phenol red and gentamicin

MATERIALS NEEDED BUT NOT PROVIDED

- » Incubator at 37°C (6% CO₂)
- » Petri dishes (e.g. Falcon IVF culture dishes (BD), Nunc culture dishes for IVF (Thermo Scientific)) and test tubes
- » Mineral oil (e.g. FertiCult™ Mineral Oil)
- » Laminar flow bench (ISO 5 environment) and microscope

PRODUCT SPECIFICATIONS

- » Chemical composition
 - » pH: 7.20-7.50 (37°C – 6% CO₂)
 - » Osmolality: 270-290 mOsm/kg
- » Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: No growth
- » Endotoxin test by Limulus Amebocyte lysate (LAL) methodology (USP <85>): < 0.25 EU/ml
- » Mouse Embryo Assay (expanded blastocysts within 96h after mouse IVF in test medium): ≥ 80%
- » Use of Ph Eur or USP grade raw materials if applicable
- » Certificate of analysis and MSDS are available upon request

WARNINGS BEFORE USE

- » Do not use the product if:
 - it becomes discoloured (if medium contains phenol red), cloudy or shows any evidence of microbial contamination
 - seal of the container is opened or defect when the product is delivered
 - expiry date has been exceeded
- » Do not freeze before use
- » Do not re-sterilize after opening
- » Products that include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.

STORAGE AND CONSERVATION

- » Store FertiCult™ IVF medium with gentamicin between 2-8°C. Store FertiCult™ IVF medium without gentamicin between 2-25°C. Once opened store all products between 2-8°C.
- » Keep away from sunlight.
- » The products can be used up to 7 days after opening, when sterile conditions are maintained and the product is stored at 2-8°C.
- » Stable after transport (max. 5 days) at elevated temperature (≤ 37°C).

Symbols as defined in ISO 15223-1:2016
Symboles tels que définis dans ISO 15223-1:2016
Symbole nach der Definition von ISO 15223-1:2016
Simboli di cui alla norma ISO 15223-1:2016
Символы соответствующие норме ISO 15223-1:2016
Simbolos conforme definido na norma ISO 15223-1:2016
Symbolen volgens ISO 15223-1:2016
Σμβολα όπως ορίζονται στο ISO 15223-1:2016

METHOD

MEDIUM PRE-EQUILIBRATION

Incubate the medium (in culture dishes or flasks with loosened screw tops) in a CO₂ incubator at 37°C and 6% CO₂. Ideally the medium should be pre-incubated overnight before use (or at least 4 hours).

Important: Read the note on optimal pH, CO₂ and atmospheric pressure below carefully.

When pre-equilibration is performed in culture dishes, follow the instructions below:

- » For micro-droplets: 6-8 droplets of 25-250 μl of FertiCult™ IVF medium may be dispensed around the 6 cm culture dish/well. To maintain pH, temperature and osmolality, make sure that the droplets are completely covered with oil (e.g. FertiCult™ Mineral Oil, FertiPro N.V.). 1-5 oocytes are usually placed per micro-droplet.
- » Alternatively, the medium can also be used for fertilization in an open system such as with the Falcon / Nunc dishes. About 1 ml of medium is placed in each well. Additional medium is placed in the reservoir surrounding the wells. This helps to maintain humidity as well as providing medium which can be used for washing of the oocytes.
- » Up to 5 oocytes per dish/well may be cultured together. In the open system, paraffin oil is not necessary, but sometimes a 1 ml layer of mineral oil may be added over the medium in the inner well. When using such an open system the incubator must be humidified.

WASHING/HOLDING OF OOCYTES AND FERTILIZATION BY IVF OR ICSI

After equilibration of the medium, the prepared dishes are ready for use. During oocyte collection, oocytes are identified in the follicular aspirates and then washed before placing into the wells or droplets in the dish. After oocyte collection, IVF oocytes will be cultured and may be rewashed, before being placed into fresh droplets/wells for further culture. Generally, within 6 hours of oocyte collection, oocytes will be inseminated by IVF or ICSI with prepared sperm and then left overnight in the incubator.

The presence of 2 pro-nuclei (2PN) is checked 16-20 hrs after fertilization. Once embryos have been identified by the presence of 2PN, they are placed into fresh dishes containing a fresh embryo culture medium (e.g. GAIN™ medium, FertiPro N.V.) for further embryo culture and embryo transfer.

NOTE: Test tube culture is now rare, but for those still using this method, place 1 ml of FertiCult™ IVF medium in each tube and then follow the same procedure as for Falcon or Nunc dishes.

NOTE ON OPTIMAL PH, CO₂ AND ATMOSPHERIC PRESSURE

For optimal embryo culture, FertiCult™ IVF medium should be used at a pH of 7.20-7.35. Different factors like the CO₂ concentration in the incubator and atmospheric pressure (which decreases at higher altitudes) have an effect on the pH after equilibration in the incubator . Therefore, we strongly advise to measure the pH under culture conditions with 6% CO₂ in the incubator in order to determine the CO₂ concentration settings of the incubator which result in an optimal pH of 7.20-7.35

PRECAUTIONS

- » Aseptic technique should be used to avoid possible contamination, even when the product contains gentamicin.
- » All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- » Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens.
- » There are no reports of proven viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- » Always wear protective clothing when handling specimens.

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