FertiCult™ IVF medium

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
LT	puslapis	11
RO	pagină	12
TR	Sayfa	13
HU	oldal	14

FertiCult™ IVF medium

STERILE A

Document ID: FP09 I07 R01 F.5
Update: 30/10/2025

MATERIAL INCLUDED

Catalogue number		
FertiCult™ IVF medium		
FECU020	5 x 20ml	FertiCult™ IVF medium
FECU050	5 x 50ml	FertiCult™ IVF medium
FECU100	3 x 100ml	FertiCult™ IVF medium
FertiCult™ IVF medium with phenol red		
FECU020PHR	5 x 20ml	FertiCult™ IVF medium with phenol red
FECU050PHR	5 x 50ml	FertiCult™ IVF medium with phenol red
FECU100PHR	3 x 100ml	FertiCult™ IVF medium with phenol red
FertiCult™ IVF medium with phenol red and gentamicin		
FECU020PHR_G	5 x 20ml	FertiCult™ IVF medium with phenol red and gentamicin
FECU050PHR_G	5 x 50ml	FertiCult™ IVF medium with phenol red and gentamicin
FECU100PHR_G	3 x 100ml	FertiCult™ IVF medium with phenol red and gentamicin

¹Conaghan J (2014). pH control in the embryo culture environment. In: Quinn P, editors. Culture media, solutions, and systems in Human ART: Cambridge University Press. P. 142-154

CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

GENERAL

Intended use

FertiCult IVF media are used for washing and holding of human oocytes, performing fertilization by in vitro fertilization (IVF), intra-cytoplasmic sperm injection (ICSI) (until 2PN stage).

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

COMPOSITION

FertiCult™ IVF medium is a ready-to-use bicarbonate-buffered balanced salt solution, supplemented with glucose (5.55mM), lactate, pyruvate, and human serum albumin (4.0g/l; medicinal substance derived from human blood plasma). The product is also available with phenol red and/or gentamicin (10mg/l; medicinal substance).

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 – 7.50 (37 °C – 6% CO₂)
- Osmolality: 270 - 290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate (LAL) methodology (USP <85>): < 0.25EU/ml
- Mouse Embryo Assay (% blastocysts after 96 hours, exposure time to test medium: 24 hours): ≥ 80%
- Human Sperm Survival Assay (% motility compared with control after 24 hours): ≥ 80%
- Use of Ph Eur or USP grade raw materials if applicable
- Certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes discoloured (if medium contains phenol red), cloudy, or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Product which include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/ warming cycles of the medium. Discard excess (unused) media.
- Keep in its original packaging until the day of use

METHOD

MEDIUM PRE-EQUILIBRATION / DISH PREPARATION

Incubate the medium (in culture dishes or flasks with loosened screw tops) in a CO₂ incubator at 37 °C and 5-6% CO₂. Ideally the medium should be pre-incubated overnight before use (or at least 5 hours), to obtain an optimal pH for oocyte washing/holding.

Important: Read the note on optimal pH, CO₂ and atmospheric pressure below carefully.

When pre-equilibration is performed in culture dishes, follow the instructions below:

- For micro-droplets: 6-8 droplets of 25-250µl of FertiCult™ IVF medium may be dispensed around the 6cm culture dish/well. To maintain pH, temperature and osmolality, make sure that the droplets are completely covered with oil (e.g. FertiCult™ Mineral Oil, FertiPro NV). 1-5 oocytes are usually placed per micro-droplet.
- Alternatively, the medium can also be used for fertilization in an open system such as with the Falcon / Nunc dishes. About 1ml of medium is placed in each well. Additional medium is placed in the reservoir surrounding the wells. This helps to maintain humidity as well as providing medium which can be used for washing of the oocytes. Up to 5 oocytes per dish/well may be cultured together. In the open system, paraffin oil is not necessary, but sometimes a 1ml layer of mineral oil may be added over the medium in the inner well. When using such an open system the incubator must be humidified.

Washing/holding of oocytes and fertilization by IVF or ICSI

After equilibration of the medium, the prepared dishes are ready for use. During oocyte collection, oocytes are identified in the follicular aspirates and then washed before placing into the wells or droplets in the dish. After oocyte collection, oocytes will be checked and may be rewashed, before being placed into fresh droplets/wells for further culture. Generally, within 6 hours of oocyte collection, oocytes will be inseminated by IVF or ICSI with prepared sperm and then left overnight in the incubator. Once embryos have been identified by the presence of 2PN(in general 16-18 hours after fertilization), they are placed into fresh dishes containing a fresh embryo culture medium (e.g GAIN™ Medium, FertiPro NV) for further embryo culture and embryo transfer.

Notes:

- Test tube culture is now rare, but for those still using this method, place 1ml of FertiCult™ IVF medium in each tube and then follow the same procedure as for Falcon or Nunc dishes.
- Each laboratory should consult its own validated procedures, optimized for its individual medical program.

Note on optimal pH, CO₂ and atmospheric pressure

For optimal embryo culture, FertiCult™ IVF medium should be used at a pH of 7.20-7.35. Different factors like the CO₂ concentration in the incubator and atmospheric pressure (which decreases at higher altitudes) have an effect on the pH after equilibration in the incubator¹. Therefore, we strongly advise to measure the pH under culture conditions with 6% CO₂ in the incubator in order to determine the CO₂ concentration settings of the incubator which result in an optimal pH of 7.20-7.35.

STORAGE/DISPOSAL INSTRUCTIONS

- Store FertiCult™ IVF medium without gentamicin between 2-25 °C. Once opened: store between 2-8 °C.
- Store FertiCult™ IVF medium with gentamicin between 2-8 °C.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from (sun)light.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for FertiCult™ IVF medium describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by downloading the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1	
	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques
	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Do not resterilize
	Contains a medicinal substance
	Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745	
	CE marking by Notified Body 2797
Symbol as defined in 21 CFR 801.15	

	US Federal law restricts this device to sale by or on the order of a physician.
---	---

FertiCult™ IVF medium

STERILE A

Document ID: FP09 I07 R01 F.5
Update: 30/10/2025

MATERIAL INCLUDED

Référence catalogue		
FertiCult™ IVF medium		
FECU020	5 x 20ml	FertiCult™ IVF medium
FECU050	5 x 50ml	FertiCult™ IVF medium
FECU100	3 x 100ml	FertiCult™ IVF medium
FertiCult™ IVF medium with phenol red		
FECU020PHR	5 x 20ml	FertiCult™ IVF medium with phenol red
FECU050PHR	5 x 50ml	FertiCult™ IVF medium with phenol red
FECU100PHR	3 x 100ml	FertiCult™ IVF medium with phenol red
FertiCult™ IVF medium with phenol red and gentamicin		
FECU020PHR_G	5 x 20ml	FertiCult™ IVF medium with phenol red and gentamicin
FECU050PHR_G	5 x 50ml	FertiCult™ IVF medium with phenol red and gentamicin
FECU100PHR_G	3 x 100ml	FertiCult™ IVF medium with phenol red and gentamicin

¹Conaghan J (2014). pH control in the embryo culture environment. In: Quinn P, editors. Culture media, solutions, and systems in Human ART: Cambridge University Press. P. 142-154

SUPPORT CLIENTS-SUPPORT TECHNIQUE

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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