# FertiVit™ Cooling

# FertiVit™ Warming

## Available languages

|    |                          |    |
|----|--------------------------|----|
| EN | <a href="#">page</a>     | 2  |
| FR | <a href="#">page</a>     | 3  |
| DE | <a href="#">Seite</a>    | 4  |
| ES | <a href="#">página</a>   | 5  |
| IT | <a href="#">pagina</a>   | 6  |
| PT | <a href="#">página</a>   | 7  |
| GR | <a href="#">σελίδα</a>   | 8  |
| NL | <a href="#">pagina</a>   | 9  |
| BG | <a href="#">страница</a> | 10 |
| TR | <a href="#">Sayfa</a>    | 11 |
| NO | <a href="#">side</a>     | 12 |
| SE | <a href="#">pagina</a>   | 13 |

# FertiVit™ Cooling

# FertiVit™ Warming

STERILE A

Document ID: FP09 I46 03 R01 B.3  
Update: 22/11/2024

## MATERIAL INCLUDED

Catalogue number

FVC\_KIT FertiVit™ Cooling kit  
FWW\_KIT FertiVit™ Warming kit

## CUSTOMER-TECHNICAL SUPPORT

FertiPro NV  
Industriepark Noord 32  
8730 Beerlem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: www.fertipro.com  
E-mail: info@fertipro.com

## IGENERAL

### Intended use

FertiVit™ Cooling / FertiVit™ Warming are a set of media for vitrification and warming of human oocytes and human embryos (zygote till blastocyst).

### Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

### Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

### Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

## MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for 4-10 vitrification / warming procedures.

FertiVit™ Cooling kit (product code: FVC\_KIT): contains one bottle each of the following media:

- Product code FPI005: 5 ml Pre-incubation medium ("PI")
- Product code FVC1001: 1 ml Cooling 1 ("C1")
- Product code FVC2001: 1 ml Cooling 2 ("C2")
- Product code FVC3001: 1 ml Cooling 3 ("C3")
- Product code FVC4001: 1 ml Cooling 4 ("C4")
- Product code FVC5001: 1.5ml Cooling 5 ("C5")

FertiVit™ Warming kit (product code: FWW\_KIT): contains one bottle of each of the following media:

- Product code FWV1005: 7.5 ml Warming 1 ("W1")
- Product code FWV2001: 1 ml Warming 2 ("W2")
- Product code FWV3001: 1 ml Warming 3 ("W3")
- Product code FWV4001: 1 ml Warming 4 ("W4")
- Product code FWV5001: 1 ml Warming 5 ("W5")
- Product code FWV6001: 1 ml Warming 6 ("W6")

The media should be used in the order displayed above (the bottles may be in a different order in the box).

## MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device, preferably closed system (FertiVit™ Cooling / FertiVit™ Warming was clinically validated with the use of Vitrisafe (Nextclinics) and High Security Vitrification device (HSV, CBS). Importantly, each laboratory should validate their own procedures when using FertiVit™ Cooling / FertiVit™ Warming in combination with other vitrification devices)

## COMPOSITION

FertiVit™ Cooling and FertiVit™ Warming are ready to use. The media are HTF-based and contain HEPES and human serum albumin (12-20 g/liter depending on the individual medium; medicinal substance derived from human blood plasma). None of the media contain antibiotics.

FertiVit™ Cooling media have increasing concentrations of DMSO and Ethylene Glycol (EG). FertiVit™ Cooling 5 medium also contains ficoll.

FertiVit™ Warming media have decreasing concentrations of sucrose.

## PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 - 7.50 (release criteria: 7.20-7.40)
- Osmolality (mOsm/kg):
  - Pre-incubation / Warming 6: 270-295 (release criteria: 270-290)
  - Warming 3: 805-865 (release criteria: 805-850)

- Warming 4: 535-565
- Warming 5: 405-435
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: no growth
- Endotoxins test by Limulus lysate methodology (USP <85>): < 0.25 EU/ml
- Mouse embryo assay (one-cell mouse embryo development until blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDs are available upon request or can be downloaded from our website (www.fertipro.com)

## WARNINGS BEFORE USE

- Do not use the product if:
  - it becomes cloudy, or shows any evidence of microbial contamination,
  - seal of the container is opened or defect when the product is delivered,
  - expiry date has been exceeded,
  - Do not freeze before use,
  - Do not re-sterilize after opening.

## STORAGE /DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

## METHOD

We strongly advice to read through all the steps of the vitrification / warming procedure before starting the procedure.

### Preliminary steps

- Mix all bottles of the kit well and warm to room temperature (20-25 °C) before use. Warming 1 might be preheated to 37 °C before use.
- Open the necessary number of vitrification devices (for maximum volume load: check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of oocytes and embryos as described in [figure 1](#).

### Notes:

- Do not use the same droplets / media set-up for different patients.
- Set-up with 6-well dish (for cases with multiple embryos/oocytes): up to 5 vitrification cycles of the same patient can be performed in the same media
- Set-up with droplets (for cases with few embryos/oocytes): the number of embryos / oocytes that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.
- Cooling 1 and 2 are only needed for vitrification of oocytes.
- Warming 5 is only needed for warming of zygotes and oocytes

### Cooling protocol

- Expose oocytes / embryos sequentially to the following media as described in [figure 2](#).


- (\*) Recommended pipetting technique in C5 ([figure 3](#)):
  - Place the oocyte / embryo on the bottom of the 1st drop of C5 (step A); Make sure the oocyte / embryo is in the drop. Preferably, do not extrude the complete volume of the C4. The oocyte / embryo will float.
  - Empty the pipette outside the drop (step B1).
  - Rinse the pipette with medium from drop 2 of C5 (step B2).
  - Aspirate a small volume of drop 2, aspirate the oocyte from drop 1, transfer to the bottom of drop 2 (step C).
  - Make sure the oocyte / embryo is in the drop. Flush in and out 2 times.
  - Place the oocyte / embryo in drop 3 and flush in and out 3 times until the shadow-like outside the drop (step D).
- Now load the oocyte / embryo on the vitrification device.

### Warming protocol


- Expose oocytes / embryos sequentially to the following media as described in [figure 4](#).

## PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.

|                     | PI      | C1    | C2    | C3    | C4    |       | C5*                                                                                                                                                                            |
|---------------------|---------|-------|-------|-------|-------|-------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Oocyte              | 1-2 min | 3 min | 3 min | 3 min | 5 - 6 | 1 min | 1 minute time to:<br>1. place the oocyte/ embryo in C5<br>2. load on the vitrification device<br>3. insert in the outer straw<br>4. seal<br>5. plunge into the liquid nitrogen |
| Zygote - 4-cell     | 1-2 min | /     | /     | 5 min | 5:30* | 1 min |                                                                                                                                                                                |
| 4-cell - blastocyst | 1-2 min | /     | /     | 5 min | 4 min | 1 min |                                                                                                                                                                                |
| DMSO / EG (%)       | 0       | 1.25  | 2.5   | 5     | 10    |       | 20                                                                                                                                                                             |

|                        | W1 (can be at 37 °C)<br>Gently swirl the straw in the medium to ensure more homogenous temperature | W2    | W3      | W4    | W5    | W6      | Pre-equilibrated culture medium |
|------------------------|----------------------------------------------------------------------------------------------------|-------|---------|-------|-------|---------|---------------------------------|
| Oocyte Zygote - 4-cell | 1 min                                                                                              | 1 min | 1-2 min | 2 min | 2 min | 1-2 min |                                 |
| 4-cell - blastocyst    | 1 min                                                                                              | 1 min | 1-2 min | 2 min | /     | 1-2 min | wash<br>cultivate               |
| sucrose (M)            | 1                                                                                                  | 0.75  | 0.5     | 0.25  | 0.125 | 0       |                                 |

- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for FertiVit™ Cooling and FertiVit™ Warming describes safety and performance characteristics for the media and is available on the website of FertiPro NV ([www.fertipro.com](http://www.fertipro.com)) or by using the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

## SYMBOLS GLOSSARY

| Symbols as defined in ISO 15223-1 |                                                                    |
|-----------------------------------|--------------------------------------------------------------------|
|                                   | Catalogue number                                                   |
|                                   | Batch code                                                         |
|                                   | Use-by date                                                        |
|                                   | Manufacturer                                                       |
|                                   | Keep away from sunlight                                            |
|                                   | Medical device                                                     |
|                                   | Consult instructions for use                                       |
|                                   | Cautions                                                           |
|                                   | Date of manufacture                                                |
|                                   | Sterilized using aseptic processing techniques                     |
|                                   | Temperature limit                                                  |
|                                   | Single sterile barrier system                                      |
|                                   | Single sterile barrier system with protective packaging outside.   |
|                                   | Do not use if package is damaged and consult instructions for use. |
|                                   | Do not resterilize                                                 |
|                                   | Contains a medicinal substance                                     |
|                                   | Contains human blood or plasma derivatives                         |

Symbol as defined in EU Commission MDR 2017/745

CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

US Federal law restricts this device to sale by or on the order of a physician.

# FertiVit™ Cooling / FertiVit™ Warming

STERILE A

Document ID: FP09 I46 03 R01 B.3  
Update: 22/11/2024

MATÉRIEL INCLUS

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|         |                       |
|---------|-----------------------|
| FVC_KIT | FertiVit™ Cooling kit |
| FWW_KIT | FertiVit™ Warming kit |

SUPPORT CLIENTS-SUPPORT TECHNIQUE

FertiPro NV  
Industriepark Noord 32  
8730 Beernem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: www.fertipro.com  
E-mail: info@fertipro.com

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- Product code FPIO05: 5 ml Pre-incubation medium ("PI")
- Product code FVC1001: 1 ml Cooling 1 ("C1")
- Product code FVC2001: 1 ml Cooling 2 ("C2")
- Product code FVC3001: 1 ml Cooling 3 ("C3")
- Product code FVC4001: 1 ml Cooling 4 ("C4")
- Product code FVC5001: 1.5ml Cooling 5 ("C5")

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- Product code FW1005: 7.5 ml Warming 1 ("W1")
- Product code FWW2001: 1 ml Warming 2 ("W2")
- Product code FWW3001: 1 ml Warming 3 ("W3")
- Product code FWW4001: 1 ml Warming 4 ("W4")
- Product code FWW5001: 1 ml Warming 5 ("W5")
- Product code FWW6001: 1 ml Warming 6 ("W6")

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