FertiVit™ Cooling

FertiVit™ Warming

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
TR	Sayfa	11
NO	side	12
SE	pagina	13
CZ	stránka	14
HU	oldal	15
DK	side	16

FertiVit™ Cooling

FertiVit™ Warming

STERILE A

Document ID: FP09 I46 03 R01 B.3
Update: 22/11/2024

MATERIAL INCLUDED

Catalogue number

FVC_KIT FertiVit™ Cooling kit
FWW_KIT FertiVit™ Warming kit

CUSTOMER-TECHNICAL SUPPORT

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

IGENERAL

Intended use

FertiVit™ Cooling / FertiVit™ Warming are a set of media for vitrification and warming of human oocytes and human embryos (zygote till blastocyst).

Indications for use

For use during ART procedures of patients and couples undergoing infertility treatments.

Intended user

The intended users are ART professionals (lab technicians, embryologists, or medical doctors).

Target patient population

The target patient population consists of patients and couples undergoing infertility treatments.

MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for 4-10 vitrification / warming procedures.

FertiVit™ Cooling kit (product code: FVC_KIT): contains one bottle each of the following media:

- Product code FPI005: 5 ml Pre-incubation medium ("PI")
- Product code FVC1001: 1 ml Cooling 1 ("C1")
- Product code FVC2001: 1 ml Cooling 2 ("C2")
- Product code FVC3001: 1 ml Cooling 3 ("C3")
- Product code FVC4001: 1 ml Cooling 4 ("C4")
- Product code FVC5001: 1.5ml Cooling 5 ("C5")

FertiVit™ Warming kit (product code: FWW_KIT): contains one bottle of each of the following media:

- Product code FW1005: 7.5 ml Warming 1 ("W1")
- Product code FW2001: 1 ml Warming 2 ("W2")
- Product code FW3001: 1 ml Warming 3 ("W3")
- Product code FW4001: 1 ml Warming 4 ("W4")
- Product code FW5001: 1 ml Warming 5 ("W5")
- Product code FW6001: 1 ml Warming 6 ("W6")

The media should be used in the order displayed above (the bottles may be in a different order in the box).

MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device, preferably closed system (FertiVit™ Cooling / FertiVit™ Warming was clinically validated with the use of Vitrifase (Nexclinics) and High Security Vitrification device (HSV, CBS). Importantly, each laboratory should validate their own procedures when using FertiVit™ Cooling / FertiVit™ Warming in combination with other vitrification devices)

COMPOSITION

FertiVit™ Cooling and FertiVit™ Warming are ready to use. The media are HTF-based and contain HEPES and human serum albumin (12-20 g/liter depending on the individual medium; medicinal substance derived from human blood plasma). None of the media contain antibiotics.

FertiVit™ Cooling media have increasing concentrations of DMSO and Ethylene Glycol (EG). FertiVit™ Cooling 5 medium also contains ficoll.

FertiVit™ Warming media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 – 7.50 (release criteria: 7.20-7.40)
- Osmolality (mOsm/kg):
 - Pre-incubation / Warming 6: 270-295 (release criteria: 270-290)
 - Warming 3: 805-865 (release criteria: 805-850)

- Warning 4: 535-565
- Warning 5: 405-435
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: no growth
- Endotoxins test by Limulus Lysate methodology (USP <85>): < 0.25 EU/ml
- Mouse embryo assay (one-cell mouse embryo development until blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination,
- seal of the container is opened or defect when the product is delivered,
- expiry date has been exceeded,
- Do not freeze before use,
- Do not re-sterilize after opening.

STORAGE / DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification / warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm to room temperature (20-25 °C) before use. Warming 1 might be preheated to 37 °C before use.
- Open the necessary number of vitrification devices (for maximum volume load; check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of oocytes and embryos as described in [figure 1](#).

Notes:

- Do not use the same droplets / media set-up for different patients.
- Set-up with 6-well dish (for cases with multiple embryos/oocytes): up to 5 vitrification cycles of the same patient can be performed in the same media
- Set-up with droplets (for cases with few embryos/oocytes): the number of embryos / oocytes that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.
- Cooling 1 and 2 are only needed for vitrification of oocytes.
- Warming 5 is only needed for warming of zygotes and oocytes

Cooling protocol

- Expose oocytes / embryos sequentially to the following media as described in [figure 2](#).

- (*) Recommended pipetting technique in C5 ([figure 3](#)):
 - Place the oocyte / embryo on the bottom of the 1st drop of C5 (step A); Make sure the oocyte / embryo is in the drop. Preferably, do not extrude the complete volume of the C4. The oocyte / embryo will float.
 - Empty the pipette outside the drop (step B1).
 - Rinse the pipette with medium from drop 2 of C5 (step B2).
 - Aspirate a small volume of drop 2, aspirate the oocyte from drop 1, transfer to the bottom of drop 2 (step C).
 - Make sure the oocyte / embryo is in the drop. Flush in and out 2 times.
 - Place the oocyte / embryo in drop 3 and flush in and out 3 times until the shadow-like shape disappears (step D).
 - Now load the oocyte / embryo on the vitrification device.

Warming protocol

- Expose oocytes / embryos sequentially to the following media as described in [figure 4](#).

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.


- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)


The SSCP for FertiVit™ Cooling and FertiVit™ Warming describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

	PI	C1	C2	C3	C4		C5*
Oocyte	1-2 min	3 min	3 min	3 min	5 - 6'	1 min	1 minute time to: 1. place the oocyte/ embryo in C5 2. load on the vitrification device 3. insert in the outer straw 4. seal 5. plunge into the liquid nitrogen
Zygote - 4-cell	1-2 min	/	/	5 min	5-30*	1 min	
4-cell - blastocyst	1-2 min	/	/	5 min	4 min	1 min	
DMSO / EG (%)	0	1.25	2.5	5	10		20

	W1 (can be at 37 °C) Gently swirl the straw in the medium to ensure more homogenous temperature	W2	W3	W4	W5	W6	Pre-equilibrated culture medium
Oocyte	1 min	1 min	1-2 min	2 min	2 min	1-2 min	
Zygote - 4-cell	1 min	1 min	1-2 min	2 min	/	1-2 min	
4-cell - blastocyst	1 min	1 min	1-2 min	2 min	/	1-2 min	
sucrose (M)	1	0.75	0.5	0.25	0.125	0	
							wash
							cultivate

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1	
	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques
	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Do not resterilize
	Contains a medicinal substance
	Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745

CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

US Federal law restricts this device to sale by or on the order of a physician.

FertiVit™ Cooling

FertiVit™ Warming

STERILE A

Document ID: FP09 I46 03 R01 B.3
Update: 22/11/2024

MATÉRIEL INCLUS

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FVC_KIT FertiVit™ Cooling kit
FWW_KIT FertiVit™ Warming kit

SUPPORT CLIENTS-SUPPORT TECHNIQUE

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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