EN FR DE IT ES

FertiVit™ Cooling / Warming kit

FertiVit™ Cooling / Warming kit

FertiVit™ Cooling/Warming kit

Media for vitrification and warming of human oocytes and embryos up to the blastocyst stage

STERILE **A**

Sterilized by sterile filtration.
Document reference: FP09 I46 03 R01 A.3
Update: 05.03.2019

INTENDED USE

FertiVit™ Cooling and FertiVit™ Warming are a set of ready-to-use media for vitrification and warming of human oocytes and embryos.

For professional use only.

BACKGROUND

Vitrification of zygotes and embryos is an established technique which has become standard practice in ART-settings. Because of European regulations defining medical safety requirements for cryopreservation of human cells; hermetically closed (aseptic) containers were developed which avoid direct contact between the embryo and liquid nitrogen during cooling and long-term storage. The FertiVit™ Cooling/Warming kit is designed to work well with the reduced cooling rates which are inherent to the use of closed vitrification devices (due to thermo-isolation). In recent years, there has also been a resurgence of interest in cryopreservation of oocytes due to the desire to preserve fertility potential in young women undergoing gonadotoxic treatments and restrictive legislation which prevent embryo cryopreservation in some countries. The American Society for Reproductive Medicine (ASRM) reported the results of a meta-analysis in order to evaluate the efficacy and safety of oocyte cryopreservation. 1200 articles were evaluated, including 4 randomized controlled trials in which the outcomes with cryopreserved and fresh oocytes in IVF/ICSI cycles were compared. From these studies, it was concluded that there is good evidence that the clinical outcome is similar for ART procedures with fresh oocytes compared to vitrified/warmed oocytes.¹
This finding has been confirmed by clinical data obtained with FertiVit™ Cooling/Warming kit.

COMPOSITION

Media are HTF-based and contain HEPES, sucrose, human serum albumin (12-20 g/liter). Cooling media also contain DMSO, Ethylene Glycol (EG) and Ficoll. None of the media contains antibiotics.

MATERIAL INCLUDED WITH THE KIT

One kit will provide sufficient medium for approximately 3-4 procedures.

FertiVit™ Cooling kit (FVC_KIT) contains one bottle each of the following media:

- » Product code FP1005: 5 ml Pre-incubation medium ("PI")
- » Product code FVC1001: 1 ml Cooling 1 ("C1")
- » Product code FVC2001: 1 ml Cooling 2 ("C2")
- » Product code FVC3001: 1 ml Cooling 3 ("C3")
- » Product code FVC4001: 1 ml Cooling 4 ("C4")
- » Product code FVC5001: 1 ml Cooling 5 ("C5")

FertiVit™ Warming kit (FWW_KIT) contains one bottle each of the following media:

- » Product code FWV1005: 5 ml Warming 1 ("W1")
- » Product code FWV2001: 1 ml Warming 2 ("W2")
- » Product code FWV3001: 1 ml Warming 3 ("W3")
- » Product code FWV4001: 1 ml Warming 4 ("W4")
- » Product code FWV5001: 1 ml Warming 5 ("W5")
- » Product code FWV6001: 1 ml Warming 6 ("W6")

The media should be used in the order displayed above (the bottles may be in a different order in the box).

MATERIAL NOT INCLUDED WITH THE KIT

- » Well dishes
- » Freezing tank with liquid nitrogen
- » Water bath (able to hold 37°C)
- » Attenuated pipettes
- » Forceps
- » Vitrification device (preferably closed device e.g. HSV straws (Cryo Bio Systems) or VitriSafe)
- » LAF-bench (ISO Class 5), Microscope, Lab timer

FERTIVIT™ COOLING/WARMING KIT AND EMBRYO CULTURE

The FertiVit™ Cooling/Warming kit can be used in combination with GAIN™, FertiCult™ IVF medium and FertiCult™ Flushing medium (FertiPro) to culture, wash oocytes and embryos before vitrification and after warming.

PRODUCT SPECIFICATIONS

- » Chemical composition
- » pH: 7,20 – 7,50 (release criteria: 7,20-7,40)
- » Osmolality (mOsm/kg):
 - » Pre-incubation / Warming 6: 270-295 (release criteria: 270-290)
 - » Warming 3: 805-865 (release criteria: 805-850)
 - » Warming 4: 535-565
 - » Warming 5: 405-435
- » Sterility: Sterile (SAL 10⁻³)
- » Endotoxins: < 0,25 EU/ml
- » Mouse Embryo Assay (blastocysts after 96h) ≥ 80%
- » Use of Ph Eur or USP grade products if applicable
- » The certificate of analysis and MSDS are available upon request

PRE-USE CHECKS

- » Do not use the product if it becomes cloudy, or shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

STORAGE INSTRUCTIONS

Store between 2-8°C. Do not freeze before use. Keep away from (sun)light. The products can be used safely up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8°C. Do not use after expiry date. Stable after transport (max. 5 days) at elevated temperature (≤37°C).

WARNINGS AND PRECAUTIONS

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid possible contamination. Only for the intended use. The long-term safety of oocyte/embryo vitrification on children born following this procedure is unknown.

METHOD

Ensure that all media bottles of the kit are well mixed before use and warmed to room temperature (20-25°C). Alternatively, Warming 1 can be warmed to 37°C. We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- » Up to 5 vitrification cycles (of the same patient) can be performed with one media set-up. Do not use the same media for different patients!
- » Open the necessary number of vitrification devices, taking into account that 1 device can hold 2-3 oocytes or 1-2 embryos, in a maximum volume load of 1µl (check the instructions for the device you are using). Conveniently place the separate parts of the device on the workbench for easy access later in the procedure.
- » Cooling procedure: In a 6-well dish, fill:
Pi: 250-300µl
C1: 250-300µl (not for embryos)
C2: 250-300µl (not for embryos)
C3: 250-300µl
C4: 250-300µl
C5: 250-300µl
- » Warming procedure: In a 6-well dish, fill:
W1: 500-800µl
W2: 250-300µl
W3: 250-300µl
W4: 250-300µl
W5: 250-300µl (not for 4-cell to blastocyst)
W6: 250-300µl

VITRIFICATION OF OOCYTES

Cooling protocol

- » Warm all media of the kit to room temperature (20-25°C) before use.
- » Oocytes are sequentially exposed to the following media:

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* Note: The complete process of placing the oocyte in "Cooling 5", loading the oocyte on the vitrification device in maximum 1µl CS, inserting the device in the outer straw and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen.

Warming protocol

- » Warm all media of the kit to room temperature (20-25°C) before use. Alternatively, Warming 1 can be warmed to 37°C (preferably in a tube).
- » In the first warming step, make sure that the straw is gently swirled in Warming medium 1 (ensures more homogenous temperature). Sequentially expose oocytes to the following media:

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* Note: Wash for 1-2 min before transfer to culture medium.

VITRIFICATION OF EMBRYOS (ZYGOTE TO BLASTOCYST)

Cooling protocol

- » Warm all media of the kit to room temperature (20-25°C) before use.
- » Embryos are sequentially exposed to the following media:

	PI	C3	C4	C5
DMSO/EG (%)	0	5	10	20
Zygotes	2 min.	5 min.	5 min. 30 sec.	40-60 sec.*
4-cell to blastocyst	2 min.	5 min.	4 min.	40-60 sec.*

* Note: The complete process of placing the embryo in "Cooling 5", loading the embryo on the vitrification device in maximum 1 µl CS, inserting the device in the outer straw and sealing should not take longer than 60 seconds before plunging the device into the liquid nitrogen.

Warming protocol

- » Warm all media of the kit to room temperature (20-25°C) before use. Alternatively, Warming 1 can be warmed to 37°C (preferably in a tube).
- » In the first warming step, make sure that the straw is gently swirled in Warming medium 1 (ensures more homogenous temperature). Sequentially expose embryos to the following media:

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