FertiCult™ Mineral Oil

FertiCult™ Mineral Oil

FertiCult™ Mineral Oil

Light mineral oil for use in IVF, ICSI and related ART procedures

STERILE A

Sterilized by sterile filtration. Document reference: FP09 I20 R01 A.17 Update: 20.12.2019

USED ABBREVIATIONS

ICSI Intracytoplasmic Sperm Injection
IVF In Vitro Fertilization
IUI Intra Uterine Insemination

GENERAL INFORMATION

FertiCult™ Mineral Oil is double washed with water, ready to use, light mineral oil (d=0.83-0.86 g/ml) for use in IVF, ICSI and related artificial reproductive techniques.

COMPOSITION

Mineral oil (paraffin oil) (EP Monograph 0240, CAS code: 8042-47-5)

INTENDED USE

Culture medium overlaid with oil was first used in the early 1960's as a technique for the culture of mammalian embryos. The oil permits small volumes of medium to be used for fertilization and embryo culture. Use of oil prevents the evaporation of media, thereby reducing the harmful effects of increased osmolality. Furthermore, the use of an oil overlay reduces temperature fluctuations and pH changes caused by loss of CO₂ from the medium when culture dishes are taken out of the incubator for embryo examination or micromanipulation. Finally, culture of multiple embryos in small volumes of medium under oil promotes the concentration of autocrine growth factors secreted by the developing embryos, leading to enhanced rates of embryo development.

For professional use only.

MATERIAL INCLUDED WITH THE KIT

| PRODUCT CODE | IN GLASS BOTTLE |
|--------------|---------------------------------|
| MINOIL050 | 1x 50ml FertiCult™ Mineral Oil |
| MINOIL100 | 1x 100ml FertiCult™ Mineral Oil |
| MINOIL500 | 1x 500ml FertiCult™ Mineral Oil |

| PRODUCT CODE | IN PETG BOTTLE |
|--------------|---------------------------------|
| MINOIL050P | 1x 50ml FertiCult™ Mineral Oil |
| MINOIL100P | 1x 100ml FertiCult™ Mineral Oil |
| MINOIL500P | 1x 500ml FertiCult™ Mineral Oil |

MATERIAL NOT INCLUDED WITH THE KIT

- » CO₂ Incubator at 37°C
- » Petri dishes
- » Cell culture media
- » LAF bench (ISO Class 5)
- » Microscope
- » Test tubes

PRODUCT SPECIFICATIONS

- » Density: 0.83-0.86 g/ml
- » Viscosity: < 30 cP at 30°C
- » Endotoxin: < 0.1 EU/ml
- » Peroxide Value (POV): < 0.1 mEq/kg
- » Sterility: SAL 10⁻³
- » Mouse embryo assay: ≥ 80% blastocysts after 120 hours in culture
- » Use of pharmaceutical grade mineral (paraffin) oil
- » A certificate of analysis and MSDS are available upon request.

PRE-USE CHECKS

- » Do not use the product if it shows any evidence of microbial contamination.
- » Do not use the product if seal of the container is opened or defect when the product is delivered.

INSTRUCTIONS FOR USE

- » FertiCult™ Mineral Oil is pre-washed. No further preparations are necessary.
- » Compatible with bicarbonate-buffered media such as GAIN medium and FertiCult™ IVF medium (FertiPro NV).
- » We advise to prepare dishes/flasks with oil the day before use. Incubate them at 37°C in a CO₂ incubator so that the oil warms up and becomes saturated with gas. If flasks are used, be sure not to tighten the cap of the flask because this will prevent gas exchange.
- » Importantly: each laboratory should consult its own optimized and validated laboratory procedures.

STORAGE AND CONSERVATION

Store FertiCult™ Mineral Oil between 15-25°C. Storage temperatures between 2-15°C are also allowed, but this might cause some turbidity which disappears when the oil is warmed up again. The product is resistant to accidental exposure to elevated temperatures during transport (≤ 37°C). Do not freeze before use. Keep the product in the original packaging or shield from light while not in use. The product can be used safely for at least 28 days after opening the container when sterile conditions are maintained and the products are stored at room temperature. Do not use after expiry date. Keep away from (sun)light and avoid exposure to sources of irradiation. Content cannot be re-sterilized after opening.

WARNINGS AND PRECAUTIONS

All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis. Always wear protective clothing when handling specimens. Always work under strict hygienic conditions (e.g. LAF-bench ISO Class 5) to avoid contamination. The oil may appear cloudy if removed from the incubator and allowed to stand at room temperature for extended periods of time. Such cloudiness usually disappears if the oil is re-equilibrated by replacing it in the incubator for several hours, if it does not disappear discard the product.

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| CODICE DEL PRODOTTO | IN FLACONE IN PETG |
|---------------------|---------------------------------|
| MINOIL050P | 1x 50ml FertiCult™ Mineral Oil |
| MINOIL100P | 1x 100ml FertiCult™ Mineral Oil |
| MINOIL500P | 1x 500ml FertiCult™ Mineral Oil |

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- Brinster R.L. (1963). A method for in vitro cultivation of mouse ova from two-cell blastocyst. *Experimental Cell Research*, 32, 205-208.
- Johnson C., Hofmann G., & Scott R. (1994). The use of oil overlay for in vitro fertilization and culture. *Assisted Reproductive Reviews*, 4, 198-201.
- Kane M.T. (1987). Culture media and culture of early embryos. *Theriogenology*, 27, 49-57.
- Swain J.E., Cabrera L., Xu X. & Smith G.D. (2012). Microdrop preparation factors influence culture-media osmolality, which can impair mouse embryo preimplantation development. *Reproductive BioMedicine Online*, 24, 142-147.
- Brinster R.L. (1966). In vitro culture of mammalian embryos. *Journal of Animal Science*, 27 (S1), 1-14.
- Jin Y., Guo X.Z., Li L., Xie C.Y., & Tan L.L. (2001). The effect of autocrine factors on development of early embryos of mouse. *Shi Yan Sheng Wu Xue Bao*, 34, 77-80.
- Provo M.B., Herr C. (1998). Washed paraffin oil becomes toxic to mouse embryos upon exposure to sunlight. *Theriogenology*, 49, 214.
- Eertmans F., Bogaert V., Volckaert D., Puype B. (2013). Validation of a potentiometric peroxide value (POV) assay for analysis of mineral oil with low oxidative content. *Journal of chemical and pharmaceutical Research* 5(11): 395-402.

TECHNICAL SUPPORT ASSISTANCE TECHNIQUE TECHNISCHE UNTERSTÜTZUNG SUPPORTO TECNICO ASSISTENCIA TÉCNICA SUPORTE TÉCNICO TECHNISCHE ONDERSTEUNING ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ

FertiPro N.V.
Industriepark Noord 32
8730 Beernem – Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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| | |
|------------|---------------------------------|
| MINOIL050P | 1x 50ml FertiCult™ Mineral Oil |
| MINOIL100P | 1x 100ml FertiCult™ Mineral Oil |
| MINOIL500P | 1x 500ml FertiCult™ Mineral Oil |

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