
10% PVP in FertiCult™ Flushing medium

STERILE A

Document ID: FP09 I17 R01 E.4
Update: 09/07/2024

MATERIAL INCLUDED/
MATERIEL INCLUS/
IM LIEFERUMFANG ENTHALTENE MATERIALIEN/
MATERIAL INCLUIDO/
MATERIALI INCLUSI/
MATERIAL INCLUIDO/
ΥΛΙΚΟ ΠΟΥ ΠΕΡΙΛΑΜΒΑΝΕΤΑΙ/
MATERIAAL AANWEZIG

CUSTOMER-TECHNICAL SUPPORT
SUPPORT CLIENTS-SUPPORT TECHNIQUE
KUNDENDIENTST-TECHNISCHER SUPPORT
ATTENCIÓN AL CLIENTE-ASISTENCIA TÉCNICA
ASSISTENZA CLIENTI-SUPPORTO TECNICO
APOYO AL CLIENTE-TÉCNICO
ΥΠΟΣΤΗΡΙΞΗ ΠΕΛΑΤΩΝ-ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ
KLANTEN-TECHNISCHE ONDERSTEUNING

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

2797

EN

10% PVP in FertiCult™ Flushing medium

INDICATIONS FOR USE

10% PVP in FertiCult™ Flushing medium is a viscous medium containing polyvinylpyrrolidone (PVP) used for sperm preparation for Intracytoplasmic Sperm Injection (ICSI). These procedures require the capture of individual sperm cells in a glass pipette for injection into the oocyte and this is facilitated by first immobilizing the sperm by placing them in a viscous medium like 10% PVP in FertiCult™ Flushing medium prior to nicking the tail to immobilize sperm completely.

The product is a ready-to-use medium.

For professional use only.

COMPOSITION

10% (w/v) dialysed Ph. Eur. grade polyvinylpyrrolidone (PVP) in FertiCult™ Flushing medium. FertiCult™ Flushing medium is an aqueous solution containing HEPES, bicarbonate, physiologic salts, glucose, lactate, pyruvate and human serum albumin (4.0g/l; substance derived from human blood plasma).

GENERAL INFORMATION

The medium is complete and needs no further additives. The medium contains HEPES: **no CO₂ incubation is required**. 10% PVP in FertiCult™ Flushing medium does not contain any antibiotics.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20-7.90 (release criteria: 7.20-7.60)
- Osmolarity: 280-340 mOsm/kg (release criteria: 280-320 mOsm/kg)
- Viscosity: >170 cP at 25 °C
- Endotoxin test by Limulus Ameboocyte Lysate Methodology (USP <85>): <0.5EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96 hours, exposure time to test medium: 10 minutes): ≥ 80%
- Human sperm survival assay (% motility compared with control, exposure time to test medium: 60 minutes): ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
 - seal of the container is opened or defect when the product is delivered.
 - expiry date has been exceeded.
 - Do not freeze before use.
 - Do not re-sterilize after opening.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Make sure that the recipient is tightly closed and avoid evaporation. Discard excess (unused) media. If multiple openings are needed, open PVPO.2S bottles via the screw cap, instead of punching the cap with a needle.
- Evaporation of the PVP solution on the slide or dish and the presence of human albumin may occasionally cause a small amount of precipitation in the medium.
- Keep in its original packaging until the day of use

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for 10% PVP in FertiCult™ Flushing medium describes safety and performance characteristics for the medium and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact our FertiPro NV for customer or technical support.

METHOD

Each laboratory should consult its own validated procedures, optimized for its individual medical program.

Standard procedure

Remove the PVP medium and preferred holding medium from storage at 2-8 °C and bring at room temperature.

1. Prepare dish:

- Depending on the number of oocytes for injection, pipette small drops of the holding medium onto the bottom of the dish (e.g. in a concentric pattern).
- Place a drop of PVP solution in the middle of the dish. The drop volume depends on the preference of the lab (e.g. 5-10µl small drop or elongated long drop up to 200µl)

Note 1:

Viscosity of the PVP solution can be diminished according to the needs in a specific ICSI procedure (e.g. when dealing with very low motility sperm) by first making a central drop with FertiCult™ Flushing medium and replacing e.g. 70% of the drop volume with PVP medium.

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