10% PVP in FertiCult™ Flushing medium

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
LT	puslapis	11
RO	pagină	12
TR	Sayfa	13
DK	side	14
HU	oldal	15

10% PVP in FertiCult™ Flushing medium

INDICATIONS FOR USE

10% PVP in FertiCult™ Flushing medium is a viscous medium containing polyvinylpyrrolidone (PVP) used for sperm preparation for Intracytoplasmic Sperm Injection (ICSI). These procedures require the capture of individual sperm cells in a glass pipette for injection into the oocyte and this is facilitated by first immobilizing the sperm by placing them in a viscous medium like 10% PVP in FertiCult™ Flushing medium prior to nicking the tail to immobilize sperm completely.

The product is a ready-to-use medium.

For professional use only.

COMPOSITION

10% (w/v) dialysed Ph. Eur. grade polyvinylpyrrolidone (PVP) in FertiCult™ Flushing medium. FertiCult™ Flushing medium is an aqueous solution containing HEPES, bicarbonate, physiologic salts, glucose, lactate, pyruvate and human serum albumin (4.0g/l, medicinal substance derived from human blood plasma).

GENERAL INFORMATION

The medium is complete and needs no further additives. The medium contains HEPES; **no CO₂ incubation is required**. 10% PVP in FertiCult™ Flushing medium does not contain any antibiotics.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20-7.90 (release criteria: 7.20-7.60)
- Osmolality: 280-340 mOsm/kg (release criteria: 280-320 mOsm/kg)
- Viscosity: >170 cP at 25 °C
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): <0.5EU/ml
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96 hours, exposure time to test medium: 10 minutes): ≥ 80%
- Human sperm survival assay (% motility compared with control, exposure time to test medium: 60 minutes): ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
 - seal of the container is opened or defect when the product is delivered.
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Make sure that the recipient is tightly closed and avoid evaporation. Discard excess (unused) media. If multiple openings are needed, open PVP0.2S bottles via the screw cap, instead of punching the cap with a needle.
- Evaporation of the PVP solution on the slide or dish and the presence of human albumin may occasionally cause a small amount of precipitation in the medium.
- Keep in its original packaging until the day of use

METHOD

Each laboratory should consult its own validated procedures, optimized for its individual medical program.

Standard procedure

Remove the PVP medium and preferred holding medium from storage at 2-8 °C and bring at room temperature.

- 1 Prepare dish:
 - a Depending on the number of oocytes for injection, pipette small drops of the holding medium onto the bottom of the dish (e.g. in a concentric pattern).
 - b Place a drop of PVP solution in the middle of the dish. The drop volume depends on the preference of the lab (e.g. 5-10µl small drop or elongated long drop up to 200µl)

Note 1:
Viscosity of the PVP solution can be diminished according to the needs in a specific ICSI procedure (e.g. when dealing with very low motility sperm) by first making a central drop with FertiCult™ Flushing medium and replacing e.g. 70% of the drop volume with PVP medium.

Note 2:
After adding PVP to the dish: immediately cover with mineral oil. Avoid evaporation of the PVP solution!
 - c Place the dish at 37 °C for 30 minutes
- 2 Add a small volume (1µl-2µl) of washed sperm into the centre of the PVP droplet.
- 3 Incubate for a (few) minute(s) to allow sperm to migrate to the periphery of the droplet.
- 4 Select the spermatozoa for injection and nick (break) the tail of the spermatozoon with the tip of the glass pipette.
- 5 Transfer the spermatozoon into one of the holding medium droplets and wash by transferring the sperm cell in and out of the medium several times.
- Suck the sperm cell into the glass pipette and use for ICSI procedure.

STORAGE/DISPOSAL INSTRUCTIONS

- Store 10% PVP in FertiCult™ Flushing medium between 2-8 °C.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from (sun)light.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic technique should be used to avoid possible contamination.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for 10% PVP in FertiCult™ Flushing medium describes safety and performance characteristics for the medium and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact our FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1

	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques

	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Do not resterilize
	Contains a medicinal substance
	Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745

CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

US Federal law restricts this device to sale by or on the order of a physician.

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