**Sil-Select Stock™** **EN**

**STERILE A**

Document ID: FP09 I13 R01 D.1  
Mise à jour: 15.11.2019

**USED ABBREVIATIONS**

EBSS Earle's Balanced Salt Solution  
ICSI Intracytoplasmic Sperm Injection  
IVF In Vitro Fertilization  
IU Intra Uterine Insemination

**INTENDED USE AND COMPOSITION**

Sil-Select Stock is a stock solution for semen preparation. It is a ready-to-use isotonic gradient for semen preparation with a density of approximately 1.2 g/ml.

Sil-Select Stock consists of saline-coated colloidal silica particles suspended in HEPEs-buffered EBSS Earle's balanced salt solution. The product code is also supplemented with 10 mg/ml gentamicin. Sil-Select Stock can be used in combination with IU, IVF and ICSI.

**MATERIAL INCLUDED**

Product code: SIS100  
 • 1x 100ml Sil-Select Stock  
 Product code: SIS500  
 • 1x 500ml Sil-Select Stock  
 Product code: SIS100\_32G  
 • 1x 100ml Sil-Select Stock with gentamicin  
 Product code: SIS500G  
 • 1x 500ml Sil-Select Stock with gentamicin

**MATERIAL NEEDED BUT NOT PROVIDED**

- 3cc syringes with 1.127 g/ml needle
- Centrifuge (must be able to operate up to 30 minutes @ 1100 g)
- Incubator at water bath at 37°C (optional)
- FertilCat® Flushing medium or Sil-Select Plus Washing/Insemination medium (FerPro® N.V.)
- LAF bench (ISO Class 5)

**PRODUCT SPECIFICATIONS**

- Chemical composition
- pH release criteria: 7.20-7.60
- pH stability criteria: 7.20-9.00
- Osmolality: 300-330 mOsm/kg
- Density: 1.150 - 1.250 g/ml
- Viscosity: < 1.75 cP
- Endotoxin level: Limulus Amoebocyte Lysate methodology (LAL) < 0.5 EU/ml
- Sterility: by the current Ph. Eur. 2.6.1. Not provided
- Spores Survived test: < 5% survival after 4 hours exposure of directly settled spermatozoa to the test medium
- Sterility by the current Ph. Eur. 2.6.1. Not provided
- A certificate of analysis and MSDS are available upon request.

**CALCULATIONS OF G-FORCES**

The g-force of your centrifuge can be calculated using this formula:

$$g = 1.118 \times r \times rpm^2 / rpm = \text{Square root } g \text{ (} \sqrt{1.118 \times r} \text{)}$$

r = radius of centrifuge in mm  
 rpm = rotations per minute / 1000

Example 1  
 r = 100 mm  
 rpm = 1800 rotations per minute  
 $g = 1.118 \times 100 \times 3.24 = 362g$

Example 2  
 r = 100 mm  
 g = 350g  
 $rpm = \sqrt{350 / 1.118} \times 1000 = 1.77 \times 1000 = 1770 \text{ rotations per minute}$

**WARNINGS BEFORE USE**

- Do not use the product if it becomes cloudy or shows any evidence of microbial contamination
- If the container is opened or defect when the product is delivered
- expiry date has been exceeded
- Do not freeze before use
- Do not re-sterilize after opening
- Products that include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics

**INSTRUCTIONS FOR PREPARATION OF GRADIENTS**

Mix the density gradient bottles by 5 bottle inversions before use.

We advise to prepare a dual gradient system (45% - 90% or 40% - 80%) starting from Sil-Select Stock. It preferred a multi-layer can be used as well (e.g. 45% - 70% - 90%). To prepare a 90% gradient, mix 1 part FertilCat® Flushing medium (or Sil-Select Plus Washing/Insemination medium) with 9 parts Sil-Select Stock.

A 45% gradient is prepared by mixing 5 parts FertilCat® Flushing medium (or Sil-Select Plus Washing/Insemination medium) with 4.5 parts of Sil-Select Stock. Alternatively any HEPEs-buffered EBSS-based medium can be used for the preparation of the gradients.

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**PRÉCAUTIONS**

- Always work under strict hygienic conditions in a LAF-bench (ISO Class 5) to avoid possible contamination. Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- All human organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.

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