EN FR DE IT ES PT NL GR

Sil-Select Plus™

STERILE A

Document ID: FP09 I13.R01.D.1
Date: 15.11.2019

Sil-Select Plus™

USED ABBREVIATIONS

ESS: Semen's balanced Salt Solution
ICSI: Intracytoplasmic Sperm Injection
IVF: In Vitro Fertilization
Intra Uterine Insemination

INTENDED USE AND COMPOSITION

Sil-Select Plus is a ready-to-use gradient system for semen preparation. The gradient system is available in a sterile disposable centrifuge tube.

Sil-Select Plus™

INSTRUCTIONS FOR USE WITH FRESH SEMEN SAMPLES

Mix the densely gradient bottles by 5 bottle inversions before use.

- Bring all components of the system and samples to room temperature or to 37°C.
- Transfer 1 ml of Sil-Select Plus Upper layer to a sterile disposable centrifuge tube.

Sil-Select Plus™

MATERIAL INCLUDES

Product code: SP008 (6 procedures)
8 x 2.5ml Sil-Select Plus UPPER LAYER
8 x 2.5ml Sil-Select Plus LOWER LAYER
8 x 2.5ml Sil-Select Plus WASHING/INSEMINATION medium

Sil-Select Plus™

INSTRUCTIONS FOR USE WITH FROZEN SEMEN SAMPLES

Mix the densely gradient bottles by 5 bottle inversions before use.

- Bring all components of the kit and samples to room temperature or to 37°C.
- Transfer 1 ml of Sil-Select Plus Upper layer to a sterile disposable centrifuge tube.

Sil-Select Plus™

MATERIAL NEEDED BUT NOT PROVIDED

Single bottles of Sil-Select Plus and Sil-Select Plus SpERM WASHING/INSEMINATION medium are all CE marked (Notified Body number 0344) and can be stored for up to 30 months, 500 µg will help to separate the sperm.

Sil-Select Plus™

PRODUCT SPECIFICATIONS

Sil-Select Plus Upper layer/Lower layer/80% media

- Chemical composition
- pH release criteria: 7.20-9.60
- Stability criteria: 2.30-9.00

Sil-Select Plus™

CALCULATIONS OF FORCES

$f = 118.1 \times r \times \rho \times \omega^2 \times r_m = \text{Squares root } f \text{ (118.1)} \times r \times \rho \times \omega^2 \times r_m$

Example 1
 $r = 100 \text{ mm}$
 $\rho = 1800 \text{ rotations per minute}$
 $f = 118.1 \times 100 \times 1.24 \times 362 = 5091$

Sil-Select Plus™

TECHNICAL SUPPORT ASSISTANCE

For further information, please contact our technical support team.

Sil-Select Plus™

FOR FURTHER INFORMATION

For further information, please contact our technical support team.

WARNINGS BEFORE USE

Do not use the product if:
- It becomes cloudy (if medium contains phenol red), cloudy or shows any evidence of microbial contamination.
- It contains visible clots or particles.
- It is past its expiration date.

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Document ID: FP09 I13 R01 D.1
Mise à jour : 15.11.2019

USED ABBREVIATIONS

EBSS Earle's Balanced Salt Solution
ICSI IntraCytoplasmic Sperm Injection
IVF In Vitro Fertilization
IU Intra Uterine Insemination

INTENDED USE AND COMPOSITION

Sil-Select Stock is a stock solution for semen preparation. It is a ready-to-use isotonic gradient for semen preparation with a density of approximately 1.2 g/ml. Sil-Select Stock consists of nitrate-coated colloidal silica particles suspended in HEPES-buffered EBSS Earle's balanced salt solution. The product code is also equipped with 10 mg/ml gentamicin. Sil-Select Stock can be used in combination with IU, IVF and ICSI.

MATERIAL INCLUDED

Product code: SIS100
 ▶ 1x 100ml Sil-Select Stock
 Product code: SIS100
 ▶ 1x 100ml Sil-Select Stock
 Product code: SIS100_326
 ▶ 1x 100ml Sil-Select Stock with gentamicin
 Product code: SIS100G
 ▶ 1x 500ml Sil-Select Stock with gentamicin

MATERIAL NEEDED BUT NOT PROVIDED

▶ 3cc syringes with 1.127 g/ml needle
 ▶ Centrifuge (must be able to operate up to 30 minutes @ 1100 g)
 ▶ Incubator at water bath at 37°C (optional)
 ▶ FerticCut Flushing medium or Sil-Select Plus Washing/Insemination medium (FerPfluo N.V.)
 ▶ LAF bench (ISO Class 5)

PRODUCT SPECIFICATIONS

▶ Chemical composition
 ▶ pH release criteria: 7.20-7.60
 ▶ pH stability criteria: 7.20-7.90
 ▶ Osmolality: 300-330 mOsm/kg
 ▶ Density: 1.150 - 1.250 g/ml
 ▶ Viscosity: ~ 1.75 cP
 ▶ Endotoxin levels: Limulus Amoebocyte Lysate methodology (LAL) < 0.5 EU/ml
 ▶ Sterility by the current Ph. Eur. 2.6.1: Not supported
 ▶ Spenn Survival test: ≥ 50% survival after 4 hours exposure of directly separated spermatozoa to the test medium.
 ▶ Sterility by the current Ph. Eur. 2.6.12: Not supported
 ▶ Certificate of analysis and MSDS are available upon request.

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

$$g = 1.118 \times r \times \text{rpm}^2 \div \text{rpm} = \text{Square root } g \div ((1.118 \times r))$$

r = radius of centrifuge in mm
 rpm = rotations per minute / 1000

Example 1
 r = 100 mm
 rpm = 1800 rotations per minute
 $g = 1.118 \times 100 \times 3.24 = 362g$

Example 2
 r = 100 mm
 g = 350g
 $\text{rpm} = \sqrt{50R / (1.118 \times 100)} = 1.77$

Example 3
 r = 100 mm
 g = 300g
 $\text{rpm} = \sqrt{450R / (1.118 \times 100)} = 1.77$

WARNINGS BEFORE USE

▶ Do not use the product if:
 - becomes cloudy or shows any evidence of microbial contamination
 - seal of the container is opened or defect when the product is delivered
 - expiry date has been exceeded
 ▶ Do not freeze before use
 ▶ Do not re-sterilize after opening
 ▶ Products that include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar aminoglycoside antibiotics

INSTRUCTIONS FOR PREPARATION OF GRADIENTS

We advise you prepare a dual gradient system (45% - 90% or 40% - 80%) starting from Sil-Select Stock. It is preferred a multi-layer can be used as well (e.g. 45% - 70% - 90%). To prepare a 90% gradient, mix 1 part FerCut Flushing medium (or Sil-Select Plus Washing/Insemination medium) with 9 parts Sil-Select Stock. A 45% gradient is prepared by mixing 5 parts FerCut Flushing medium (or Sil-Select Plus Washing/Insemination medium) with 4.5 parts of Sil-Select Stock. Alternatively, any HEPES-buffered EBSS-based medium can be used for the preparation of the gradient.

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