SpermFreeze™ SSP

INTENDED USE

SpermFreeze™ SSP is a medium for cryopreservation of human sperm for further use in Assisted Reproductive Technologies.

For professional use only.

COMPOSITION

SpermFreeze™ SSP is a ready-to-use HEPES buffered cryopreservation medium which also contains physiologic salts, glycine, glucose, lactate and the cryoprotectants glycerol (27%), sucrose and Human Serum Albumin (4.0g/ Liter; medicinal substance derived from human blood plasma) to protect the sperm from damage during the freezing procedure.

Note: the concentration of glycerol in SpermFreeze™ SSP is significantly higher than in SpermFreeze™, therefore it is important to use the correct dilution factors when using SpermFreeze™ versus SpermFreeze™ SSP.

MATERIAL NOT INCLUDED BUT NEEDED

- Sperm freezing straws (e.g. CBS™ high security spermstraws) or cryo tubes
- Freezing tank with liquid nitrogen

PRODUCT SPECIFICATIONS

- Chemical composition
- pH between 7.20 – 7.90 (Release criteria: 7.20 – 7.60)
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>: No growth
- Endotoxins (USP <85>): < 0.25EU/ml
- Human sperm survival assay (% motility compared with control after 4 hours): ≥80%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
- Keep in its original packaging until the day of use

STORAGE/DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Stable after transport (max. 5 days) at elevatedtemperature (≤37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

Ensure all media are well mixed and at room temperature before use.

Before freezing

Ideally, cryopreservation is performed on native semen samples. In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted with FertiCult™ Flushing medium before freezing, if preferred.

- Allow the semen to liquify at room temperature for 30 minutes.
- Mix 3 parts sperm with 1 part of SpermFreeze™ SSP. Importantly, add the medium in drops while gently swirling. Caution: to avoid cold-shock, make sure SpermFreeze™ SSP is at room temperature.
- Leave the mixture for 10 minutes at room temperature for equilibration.
- Suck the sample/medium mixture into the freezing straws, leaving approximately 1.5cm of air at the end of the straw.
- Make sure that the outer surface of the straw is clean (no sperm) and seal the straw.
- If applicable: dry off individually with a dust free cloth.
- Shake to move the air-bubble to the center of the straw
- Place the straws horizontally (e.g. on a styofam board) in a liquid nitrogen bath to allow for freezing in vapor phase. Leave for (at least) 15 minutes.
- Transfer straws quickly into liquid nitrogen and store at -196 °C.

Thawing

- Remove as many straws as required from the liquid nitrogen.
- Place the straws in tap water for 5 minutes (room

temperature or 37 °C).

- Cut off the end of the straw, place the open end insidea container (e.g. a test tube) and tap the straw against the side of the container to allow complete evacuation of the mixture.
- Dilute the concentrated sperm in a suitable insemination medium (e.g. FertiCult™ Flushing medium) (at least 3ml per 0.5ml semen) and mix thoroughly.
- Centrifuge during 15 minutes at 300-350g.
- Resuspend pellet in a suitable insemination medium (e.g. FertiCult™ Flushing medium) and access recovery.

Notes:


- If necessary, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium (e.g. FertiCult™ Flushing medium).
- For additional details on the method, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for its individual medical program.

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was negative for Anti-HIV-1/2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination. SpermFreeze™ SSP does not containany antibiotics.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the competent authority of the EU Member State in which the user and/ or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

Although damaging effects on sperm motility, vitality, DNA fragmentation and fertilizing potential may be inherent to the cryopreservation process, these remain within acceptable limits for SpermFreeze™ SSP as reported in scientific literature. The SSCP for SpermFreeze™ SSP describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

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