# SpermFreeze™ SSP

EN

## SpermFreeze™ SSP

### INTENDED USE

SpermFreeze™ SSP is a medium for cryopreservation of human sperm for further use in Assisted Reproductive Technologies.

For professional use only.

### COMPOSITION

SpermFreeze™ SSP is a ready-to-use HEPES buffered cryopreservation medium which also contains physiologic salts, glycine, glucose, lactate and the cryoprotectors glycerol (27%), sucrose and Human Serum Albumin (4.0g/l; medicinal substance derived from human blood plasma) to protect the sperm from damage during the freezing procedure.

**Note:** the concentration of glycerol in SpermFreeze™ SSP is significantly higher than in SpermFreeze™, therefore it is important to use the correct dilution factors when using SpermFreeze™ versus SpermFreeze™ SSP.

### ACCESSORIES

The following device can be used in combination with SpermFreeze™ SSP to further standardize the cryopreservation process:

- SpermFreeze Box, product code SPF\_BOX (FertiPro NV)

### MATERIAL NOT INCLUDED BUT NEEDED

- Sperm freezing straws (e.g. CBS™ high security spermstraws) or cryo tubes
- Freezing tank with liquid nitrogen

### PRODUCT SPECIFICATIONS

- Chemical composition
- pH between 7.20 – 7.90 (Release criteria: 7.20 – 7.60)
- Sterility test by the current Ph. Eur. 2.6.1 / USP <71>; No growth
- Endotoxins (USP <85>): < 0.25EU/ml
- Human sperm survival assay (% motility compared with control after 4 hours): ≥ 80%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request or can be downloaded from our website ([www.fertipro.com](http://www.fertipro.com)).

### WARNINGS BEFORE USE

- Do not use the product if:
  - it becomes cloudy, or shows any evidence of microbial contamination;
  - seal of the container is opened or defect when the product is delivered;
  - expiry date has been exceeded;
  - Do not freeze before use.
  - Do not re-sterilize after opening.
  - Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
  - Keep in its original packaging until the day of use

### STORAGE/DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

### METHOD

Ensure all media are well mixed and at room temperature before use.

### Before freezing

Ideally, cryopreservation is performed on native semen samples. In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted with FertiCult™ Flushing medium before freezing, if preferred.

### Freezing

- 1 Allow the semen to liquefy at room temperature for 30 minutes.

- 2 Mix 3 parts sperm with 1 part of SpermFreeze™ SSP. Importantly: add the medium in drops while gently swirling. Caution: to avoid cold-shock, make sure SpermFreeze™ SSP is at room temperature.

- 3 Leave the mixture for 10 minutes at room temperature for equilibration.

- 4 Suck the sample/medium mixture into the freezing straw, leaving approximately 1.5cm of air at the end of the straw.

- 5 Make sure that the outer surface of the straw is clean (no sperm) and seal the straw.

- 6 If applicable: dry off individually with a dust free cloth.

- 7 Shake to move the air-bubble to the center of the straw.

- 8 Place the straws horizontally (e.g. on a styrofoam board) in a liquid nitrogen bath to allow for freezing in vapor phase. Leave for (at least) 15 minutes.

- 9 Transfer straws quickly into liquid nitrogen and store at -196 °C.

When using the SpermFreeze Box as accessory, step 8 and 9 in the freezing process are adapted as follows:

- Fill the box with 1-2 liter of liquid nitrogen

- Place your sealed straws (maximum 15) on the float

- Deposit the float on the liquid nitrogen

- Close the box by placing the lid loosely on the box

- Leave for (at least) 15 minutes

STERILE A

Sterilized by aseptic filtration

Document ID: FP09 i11 SSP R01 C.4

Update: 14/06/2024

**MATERIAL INCLUDED/**  
MATERIAL INCLUS/ **LIEFERUMFANG ENTHALTNE MATERIALIEN/**  
MATERIAL INCLUIDO/ **MATERIAL INCLUIDO/**  
MATERIAL INCLUIDO/ **ΥΛΙΚΟ ΠΟΥ ΠΕΡΙΑΜΒΑΝΕΤΑΙ/**  
MATERIAL AANWEZIG

Catalogue number	Reference catalogue	Bestellnummer	Número de catálogo	Auftrags-Nr.	Catalogusnummer
SSP125	1 x 125ml	SpermFreeze™ SSP			
SSP020	1 x 20ml	SpermFreeze™ SSP			
SSP001	5 x 1ml	SpermFreeze™ SSP			

**CUSTOMER-TECHNICAL SUPPORT**  
SUPPORT CLIENTS-SUPPORT TECHNIQUE  
KUNDENDIENT-TECHNISCHE SUPPORT  
ATENCION AL CLIENTE-ASISTENCIA TECNICA  
ASSISTENZA CLIENTI-SUPPORTO TECNICO  
APOYO AL CLIENTE-TECNICO  
ΥΠΟΣΤΗΡΙΞΗ ΠΕΛΑΤΩΝ-ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ  
KLANTEN-TECHNISCHE ONDERSTEUNING

FertiPro NV  
Industrieweg Noord 32  
8730 Beernem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: [www.fertipro.com](http://www.fertipro.com)  
E-mail: [info@fertipro.com](mailto:info@fertipro.com)

2797

• Plunge the straws in the liquid nitrogen  
• Store them in your nitrogen storage container for long term preservation

### Thawing

- 1 Remove as many straws as required from the liquid nitrogen
- 2 Place the straws in tap water for 5 minutes (room temperature or 37 °C)
- 3 Cut off the end of the straw, place the open end of the straw against the side of the container to allow complete evacuation of the mixture
- 4 Dilute the concentrated sperm in a suitable insemination medium (e.g. FertiCult™ Flushing medium) (at least 3ml per 0.5ml semen) and mix thoroughly
- 5 Centrifuge during 15 minutes at 300-350g
- 6 Resuspend pellet in a suitable insemination medium (e.g. FertiCult™ Flushing medium) and access recovery.

### Notes:

- 1 If necessary, use sperm preparation techniques after thawing the semen to eliminate dead sperm cells and debris. Dilute the concentrated sperm in a suitable insemination medium (e.g. FertiCult™ Flushing medium).
- 2 For additional details on the method, each laboratory should consult its own laboratory procedures and protocols which have been specifically developed and optimized for its individual medical program.

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### WARNINGS BEFORE USE

- Do not use the product if:
  - it becomes cloudy, or shows any evidence of microbial contamination;
  - seal of the container is opened or defect when the product is delivered;
  - expiry date has been exceeded;
  - Do not freeze before use.
  - Do not re-sterilize after opening.
  - Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
  - Keep in its original packaging until the day of use

### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

Although damaging effects on sperm motility, DNA fragmentation and fertilizing potential upon thawing are inherent to the cryopreservation process, these remain within acceptable limits for SpermFreeze™ SSP. As reported in scientific literature, the cryopreservation is effected upon the selection of donors and plasma pools for specific markers of selection and the inclusion of effective freezing steps for the inactivation/removal of viruses. Despite the fact that medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infectious agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Always wear protective clothing when handling specimens.

Aseptic techniques should be used to avoid possible contamination. SpermFreeze™ SSP does not contain any antibiotics.

Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

For further questions regarding the safety and performance, please contact FertiPro NV for customer or technical support.

### STORAGE/DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

### METHOD

Ensure all media are well mixed and at room temperature before use.

### Before freezing

Ideally, cryopreservation is performed on native semen samples. In case of very low sperm concentrations it is advisable to concentrate the sperm before freezing. In case of very high sperm concentration, sperm can be diluted with FertiCult™ Flushing medium before freezing, if preferred.

### Freezing

- 1 Allow the semen to liquefy at room temperature for 30 minutes.

- 2 Mix 3 parts sperm with 1 part of SpermFreeze™ SSP. Importantly: add the medium in drops while gently swirling. Caution: to avoid cold-shock, make sure SpermFreeze™ SSP is at room temperature.

- 3 Leave the mixture for 10 minutes at room temperature for equilibration.

- 4 Suck the sample/medium mixture into the freezing straw, leaving approximately 1.5cm of air at the end of the straw.

- 5 Make sure that the outer surface of the straw is clean (no sperm) and seal the straw.

- 6 If applicable: dry off individually with a dust free cloth.

- 7 Shake to move the air-bubble to the center of the straw.

- 8 Place the straws horizontally (e.g. on a styrofoam board) in a liquid nitrogen bath to allow for freezing in vapor phase. Leave for (at least) 15 minutes.

- 9 Transfer straws quickly into liquid nitrogen and store at -196 °C.

When using the SpermFreeze Box as accessory, step 8 and 9 in the freezing process are adapted as follows:

- Fill the box with 1-2 liter of liquid nitrogen

- Place your sealed straws (maximum 15) on the float

- Deposit the float on the liquid nitrogen

- Close the box by placing the lid loosely on the box

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