VitriFreeze™

VitriThaw™

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
RO	pagină	11
TR	Sayfa	12
DK	side	13
HU	oldal	14

VitriFreeze™

VitriThaw™

Document ID: FP09 I46 R01 C.2
Update: 21/11/2023

MATERIAL INCLUDED

Catalogue number	
VF_KIT1	VitriFreeze
VT_KIT1	VitriThaw

CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

INTENDED USE

VitriFreeze™ and VitriThaw™ are a set of media for vitrification and thawing of human embryos (morula till expanded blastocyst stage).

For professional use only.

MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for a maximum of 6 vitrification / warming procedures, depending on the protocol used.

- VitriFreeze™ kit (product code VF_KIT1):
contains one bottle of the following media:
- Product code VPI005: 5 ml VitriFreeze Pre-incubation (“PI”)
 - Product code VF1001: 1 ml VitriFreeze Freezing 1 (“F1”)
 - Product code VF2001: 1 ml VitriFreeze Freezing 2 (“F2”)

- VitriThaw™ kit (product code: VT_KIT1):
contains one bottle of the following media:
- Product code VT1005: 5 ml VitriThaw Thawing 1 (“T1”)
 - Product code VT2001: 1 ml VitriThaw Thawing 2 (“T2”)
 - Product code VT3001: 1 ml VitriThaw Thawing 3 (“T3”)
 - Product code VT4001: 1 ml VitriThaw Thawing 4 (“T4”)

The media should be used in the order displayed above (the bottles may be in a different order in the kit).

MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device (Although VitriFreeze™/VitriThaw™ can be used with open and closed devices, the use of closed devices is recommended. VitriFreeze™/VitriThaw™ were clinically validated with use of the following closed devices: High Security Vitrification straw (HSV; CBS) or Vitrisafe (Nextclinics). Importantly, each laboratory should validate their own procedures when using VitriFreeze™/VitriThaw™ in combination with other vitrification devices)

COMPOSITION

VitriFreeze™ and VitriThaw™ are ready to use. The media contain PBS and human serum albumin (12-20 g/l depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

VitriFreeze™ media have increasing concentrations of DMSO and ethylene glycol (EG). VitriFreeze Freezing 2 also contains sucrose and Ficoll. VitriThaw™ media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 – 7.40
- Osmolality (mOsm/kg):
 - VitriFreeze Pre-incubation: 270-295 (release criteria: 270-290)
 - VitriThaw Thawing 1: 805-865 (release criteria: 805-850)
 - VitriThaw Thawing 2: 535-565
 - VitriThaw Thawing 3: 405-435
 - VitriThaw Thawing 4: 270-295 (release criteria: 270-290)
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: no growth
- Endotoxins test by Limulus Lysate methodology (USP <85>): < 0.25 EU/ml
- Mouse embryo assay (one-cell mouse embryo development until blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80 %
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
 - seal of the container is opened or defect when the product is delivered.
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.

STORAGE /DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification/ warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm all media to an equal temperature (room temperature (RT, 20-25 °C) or 37 °C) depending on the embryo stage or medium (see figure 2 and 4).
- Open the necessary number of vitrification devices, taking into account that 1 device can hold 1-2 embryos, (for maximum volume load: check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of embryos as described in figure 1.

Notes:

- Do not use the same droplets/ media set-up for different patients.
- Set-up with 4-well dish (for cases with multiple embryos): up to 5 vitrification cycles of the same patient can be performed in the same media.
- Set-up with droplets (for cases with few embryos): the number of embryos that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.

Vitrification protocol

- Expose the embryos sequentially to the following media as described in figure 2.

- (!) Recommended pipetting technique in F2 (figure 3):
 - Place the embryo on the bottom of the 1st drop of F2 (step A); Make sure the embryo is in the drop. Preferably, do not extrude the complete volume of the F1. The embryo will float.
 - Empty the pipette outside the drop (step B1).
 - Rinse the pipette with medium from drop 2 of F2 (step B2).
 - Aspirate a small volume of drop 2, aspirate the embryo from drop 1, transfer to the bottom of drop 2 (step C).
 - Make sure the embryo is in the drop. Flush in and out 2 times.
 - Place the embryo in drop 3 and flush in and out 3 times until the shadow-like shape disappears (step D).
 - Now load the embryo on the vitrification straw.

Thawing protocol

- Expose the embryos sequentially to the following media as defined in figure 4.

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.


Figure 2

stage		PI	F1	F2* (see figure 3)	
morulae early Blastocyst	RT	2 min	2 min	30 - 60 sec	*60 seconds time to: 1. place the embryo in F2 2. load on the vitrification device 3. insert in the outer straw 4. seal 5. plunge into the liquid nitrogen
blastocyst expanded Blastocyst	RT or 37 °C	2 min	2min (when F1 is used at RT) 3min (when F1 is used at 37 °C)	30 - 60 sec	
DMSO / EG (%)	-	0	10	20	


Figure 4

	temperature	Preferable 37 °C (or RT)	RT			Pre-equilibrated culture medium	
		T1 Gently swirl the straw in the medium to ensure more homogenous temperature	T2	T3	T4		
All embryo stages		3 min	2 min	2 min	1 min	wash	cultivate
sucrose (M)		0.5	0.25	0.125	0		

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for VitriFreeze™ and VitriThaw™ describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1	
REF	Catalogue number
LOT	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
MD	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
STERILE A	Sterilized using aseptic processing techniques
°C	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Do not resterilize
	Contains a medicinal substance
	Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745

CE 2797 CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

R ONLY US Federal law restricts this device to sale by or on the order of a physician.

VitriFreeze™

VitriThaw™

Document ID: FP09 I46 R01 C.2
Update: 21/11/2023

MATÉRIEL INCLUS

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VF_KIT1	VitriFreeze
VT_KIT1	VitriThaw

SUPPORT CLIENTS-SUPPORT TECHNIQUE

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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