FertiCult™ High Viscosity Oil

EN

FertiCult™ High Viscosity Oil

INDICATIONS FOR USE

FertiCult High Viscosity Oil is pre-washed with water, ready to use, high viscosity mineral oil ($d=0.86-0.88\text{g/ml}$) for use in In Vitro Fertilization (IVF), Intracytoplasmatic Sperm Injection (ICSI) and Assisted reproductive technologies (ART).

Oil overlay is highly recommended during oocyte culture, denudation of oocytes, fertilization of oocytes (either by IVF or ICSI), embryo culture and manipulation. The oil permits small volumes of medium to be used for fertilization and embryo culture / handling. The primary function of an oil overlay is to prevent the culture medium from evaporation, which would result in increased osmotic pressure, and to avoid fluctuations in the pH and temperature. This fact is critical since variations in the pH of the medium may affect embryo metabolism, whereas osmotic stress may cause DNA and protein alterations that affect cell function.

For professional use only.

COMPOSITION

Mineral oil (Paraffin oil) (EP Monograph 0239, CAS code: 8042-47-5)

PRODUCT SPECIFICATIONS

- Chemical composition
- Density: $0.86-0.88\text{g/ml}$
- Viscosity: $95-140\text{cP}$ at 30°C
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): $< 0.1\text{EU}/\text{ml}$
- Peroxide Value (POV): $< 0.1\text{mEq/kg}$
- Sterility test according to current Ph. Eur. 2.6.1/ USP <71>: no growth
- Mouse embryo assay: $\geq 80\%$ blastocysts after 120 hours culture in standard mouse IVF medium covered with test item
- Use of pharmaceutical grade mineral (paraffin) oil
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use.

STORAGE/DISPOSAL INSTRUCTIONS

- Store between $15-25^\circ\text{C}$.
- Keep away from (sun)light and avoid exposure to sources of irradiation.
- The products can be used up to 1 month after opening, when sterile conditions are maintained and the products are stored at room temperature.
- Stable after transport (max. 5 days) at elevated temperature ($\leq 37^\circ\text{C}$).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

- For professional use only.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- The oil may appear cloudy if removed from the incubator and allowed to stand at room temperature for extended periods of time. Such cloudiness usually disappears if the oil is re-equilibrated by replacing it in the incubator for several hours, if it does not disappear discard the product.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

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HVOIL050P	1x50ml FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL100P	1x100ml FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL500P	1x500ml FertiCult™ High Viscosity Oil (PETG bottle)
HVOIL50	1x50ml FertiCult™ High Viscosity Oil (glass bottle)
HVOIL100	1x100ml FertiCult™ High Viscosity Oil (glass bottle)
HVOIL500	1x500ml FertiCult™ High Viscosity Oil (glass bottle)

CUSTOMER-TECHNICAL SUPPORT
SUPPORT CLIENTS-SUPPORT TECHNIQUE
KUNDENDRIST-TECHNISCHER SUPPORT
ATENCION AL CLIENTE-ASISTENCIA TECNICA
ASSISTENZA CLIENTI-SUPPORTO TECNICO
APOYO AL CLIENTE-TÉCNICO
ΥΠΟΣΤΗΡΙΞΗ ΠΕΛΑΤΩΝ-ΤΕΧΝΙΚΗ ΥΠΟΣΤΗΡΙΞΗ
KLANTEN-TECHNISCHE ONDERSTEUNING

FertiPro NV
Industriepark Noord 32
8730 Beerne / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

CE
2797

PRECAUTIONS

—

- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- The oil may appear cloudy if removed from the incubator and allowed to stand at room temperature for extended periods of time. Such cloudiness usually disappears if the oil is re-equilibrated by replacing it in the incubator for several hours, if it does not disappear discard the product.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

CLINICAL BENEFIT

The clinical benefit of FertiCult High Viscosity Oil is related to the fact that the oil is used during ART-procedures of patients and couples with infertility problems. It is expected that the outcomes of ART procedures in which the oil is used are comparable with the ART outcomes (i.e. pregnancy rate and delivery rate) as published in the annual peer-reviewed report of the European IVF-monitoring Consortium (EIM) for the European Society of Human Reproduction and Embryology (ESHRE).

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

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