# FertiCult™

## Aspiration medium

OUTDATED

### Available languages

EN	<a href="#">page</a>	2
FR	<a href="#">page</a>	3
DE	<a href="#">Seite</a>	4
ES	<a href="#">página</a>	5
IT	<a href="#">pagina</a>	6
PT	<a href="#">página</a>	7
GR	<a href="#">σελίδα</a>	8
NL	<a href="#">pagina</a>	9
BG	<a href="#">страница</a>	10
LT	<a href="#">puslapis</a>	11
RO	<a href="#">pagină</a>	12
SK	<a href="#">stran</a>	13
TR	<a href="#">Sayfa</a>	14
DK	<a href="#">side</a>	15
HU	<a href="#">oldal</a>	16

# FertiCult™ Aspiration medium

**STERILE A**

Document ID: FP09 I19 R01 D.1  
Update: 22/06/2023

## MATERIAL INCLUDED

Catalogue number

ASPI020	5x 20ml FertiCult™ Aspiration medium
ASPI100	5x 100ml FertiCult™ Aspiration medium

## CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV  
Industriepark Noord 32  
8730 Beernem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: [www.fertipro.com](http://www.fertipro.com)  
E-mail: [info@fertipro.com](mailto:info@fertipro.com)

## INDICATIONS FOR USE

FertiCult™ Aspiration medium is a ready-to-use formulation used for oocyte aspiration and follicular flushing.

*For professional use only.*

## GENERAL INFORMATION

FertiCult™ Aspiration medium is complete and needs no further additives. **The medium contains HEPES and does not require CO<sub>2</sub> incubation.**

FertiCult™ Aspiration medium has to be pre-warmed to 37 °C in the incubator with lid closed (approximately 12 hours, but dependent on product volume).

## COMPOSITION

FertiCult™ Aspiration medium is an aqueous solution containing HEPES, bicarbonate, physiologic salts, glucose, lactate, pyruvate and 2.5 IU/ml heparin (medicinal substance derived from porcine origin).

## PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.30-7.90 (release criteria: 7.30-7.60)
- Osmolality: 270-290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.25 EU/ml
- One-cell mouse embryo assay (% blastocysts after 96h, exposure time to test medium: 60 minutes): ≥ 80%
- Use of Ph Eur or USP grade raw materials if applicable
- Certificate of analysis and MSDS are available upon request or can be downloaded from our website ([www.fertipro.com](http://www.fertipro.com)).

## WARNING BEFORE USE

- Do not use the product if:
  - it becomes cloudy, or shows any evidence of microbial contamination.
  - seal of the container is opened or defect when the product is delivered.
  - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple opening/warming cycles of the medium. Discard excess (unused) media.

## STORAGE/ DISPOSAL INSTRUCTIONS

- Store FertiCult™ Aspiration medium between 2-25 °C. Once opened: store between 2-8 °C.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep in its original packaging until day of use. Keep away from (sun)light.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

## METHOD

Each laboratory should consult its own validated procedures, optimized for its individual medical program.

## SUGGESTED PROTOCOL FOR OOCYTE COLLECTION

At oocyte collection, keep FertiCult™ Aspiration medium at 37 °C or as close to this as possible. Take about 20-50ml medium at a time, using a rubber free syringe (e.g. Braun Injekt). Assemble needle and tubing set according to manufacturer's instructions. Attach syringe containing FertiCult™ Aspiration medium to the appropriate luer, allowing connection to the lumen of the needle. The person holding the syringe should try to apply warmth to it by holding it in the palm of the hand covering as much of the syringes surface as possible. The tubing and needle can be flushed with FertiCult™ Aspiration medium before start of the procedure.

During aspiration of follicles, if an ovum is not extracted after a follicle has been emptied, flushing is indicated. The collapsed follicle may be refilled with FertiCult™ Aspiration medium, after which it may be aspirated once again. This may be repeated (usually three attempts), after refilling the syringe with fresh warm aspiration medium each time, until an ovum is recovered or the surgeon decides to abandon the follicle.

After oocyte collection, transfer oocytes to an appropriate culture medium (for example GAIN Medium, FertiPro NV).

## PRECAUTIONS

- FertiCult™ Aspiration medium contains heparin which is derived from porcine intestinal mucosa. The heparin is certified with a Certificate of Suitability (CEP). The animals from which the heparin is derived, are determined "fit for human consumption".
- Aseptic techniques should be used to avoid possible contamination.
- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)



The SSCP for FertiCult™ Aspiration medium describes safety and performance characteristics for the medium and is available on the website of FertiPro NV ([www.fertipro.com](http://www.fertipro.com)) or by downloading the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

## SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1

	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques

	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Contains biological material of animal origin

Symbol as defined in EU Commission MDR 2017/745

	CE marking by Notified Body 2797
---	----------------------------------

Symbol as defined in 21 CFR 801.15

	US Federal law restricts this device to sale by or on the order of a physician.
---	---

# FertiCult™

## Aspiration medium

**STERILE A**

Document ID: FP09 I19 R01 D.1  
Update: 22/06/2023

### MATÉRIEL INCLUS

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ASPI020	5x 20ml	FertiCult™ Aspiration medium
ASPI100	5x 100ml	FertiCult™ Aspiration medium

### SUPPORT CLIENTS-SUPPORT TECHNIQUE

FertiPro NV  
Industriepark Noord 32  
8730 Beernem / Belgium  
Tel +32 (0)50 79 18 05  
Fax +32 (0)50 79 17 99  
URL: [www.fertipro.com](http://www.fertipro.com)  
E-mail: [info@fertipro.com](mailto:info@fertipro.com)

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