
FertiCult™ Mineral Oil


OUTDATED

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
LT	puslapis	11
RO	pagină	12
SK	stran	13
TR	Sayfa	14
HU	oldal	15
DK	side	16
CZ	stránka	17

FertiCult™ Mineral Oil

STERILE A

Document ID: FP09 I20 R01 B.3
Update: 08/05/2023

MATERIAL INCLUDED

Catalogue number		
MINOIL050	1 x 50ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL100	1 x 100ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL500	1 x 500ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL050P	1 x 50ml	FertiCult™ Mineral Oil (PETG bottle)
MINOIL100P	1 x 100ml	FertiCult™ Mineral Oil (PETG bottle)
MINOIL500P	1 x 500ml	FertiCult™ Mineral Oil (PETG bottle)

CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

INDICATIONS FOR USE

FertiCult™ Mineral Oil is double washed with water, ready to use, light mineral oil (d=0.83-0.86g/ml) for use in In Vitro Fertilization (IVF), Intracytoplasmatic Sperm Injection (ICSI) and related Assisted reproductive technologies (ART).

Oil overlay is highly recommended during oocyte culture, denudation of oocytes, fertilization of oocytes (either by IVF or ICSI), embryo culture and manipulation. The oil permits small volumes of medium to be used for fertilization and embryo culture / handling. The primary function of an oil overlay is to prevent the culture medium from evaporation, which would result in increased osmotic pressure, and to avoid fluctuations in the pH and temperature. This fact is critical since variations in the pH of the medium may affect embryo metabolism, whereas osmotic stress may cause DNA and protein alterations that affect cell function.

For professional use only.

COMPOSITION

Mineral oil (Paraffin oil) (CAS code: 8042-47-5)

PRODUCT SPECIFICATIONS

- Chemical composition
- Density: 0.83-0.86g/ml
- Viscosity: < 60cP at 30 °C
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.1EU/ml
- Peroxide Value (POV): <0.1mEq/kg (release criteria: <0.06 mEq/kg)
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Mouse embryo assay: ≥ 80% blastocysts after 120 hours in standard mouse IVF medium covered with test item
- Use of pharmaceutical grade mineral (paraffin) oil
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use

STORAGE/DISPOSAL INSTRUCTIONS

- Store between 15-25 °C. If this storage temperature can not be ensured, FertiCult Mineral Oil can alternatively be stored at 2-8 °C. Note that this might cause some turbidity which disappears when the oil is warmed up again.
- Keep away from (sun)light and avoid exposure to sources of irradiation
- The products can be used up to 28 days after opening, when sterile conditions are maintained and the products are stored at room temperature, or alternatively at 2-8 °C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

INSTRUCTIONS FOR USE

- FertiCult™ Mineral Oil is pre-washed. No further preparations are necessary.
- Compatible with gamete/embryo culture or handling media classified as class IIb or III Medical Device (for example but not limited to: FertiCult™ IVF medium, GAIN™ medium and other products of product range of ART media manufactured by FertiPro NV)
- We advise to prepare dishes/flasks with oil the day before use. Incubate them at 37 °C in a CO₂ incubator so that the oil warms up and becomes saturated with gas. If flasks are used, be sure not to tighten the cap of the flask because this will prevent gas exchange.
- Importantly: each laboratory should consult its own optimized and validated laboratory procedures.

PRECAUTIONS

- All human, organic material should be considered potentially infectious. Handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens
- Aseptic techniques should be used to avoid possible contamination.
- The oil may appear cloudy if removed from the incubator and allowed to stand at room temperature for extended periods of time. Such cloudiness usually disappears if the oil is re-equilibrated by replacing it in the incubator for several hours, if it does not disappear discard the product.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, the competent authority of the EU Member State in which the user and/or patient is established.

CLINICAL BENEFIT

The clinical benefit of FertiCult Mineral Oil is related to the fact that the oil is used during ART-procedures of patients and couples with infertility problems. It is expected that the outcomes of ART procedures in which the oil is used are comparable with the ART outcomes (i.e. pregnancy rate and delivery rate) as published in the annual peer-reviewed report of the European IVF-monitoring Consortium (EIM) for the European Society of Human Reproduction and Embryology (ESHRE).

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1	
	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques
	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
Symbol as defined in EU Commission MDR 2017/745	
	CE marking by Notified Body 2797
Symbol as defined in 21 CFR 801.15	
	US Federal law restricts this device to sale by or on the order of a physician.

FertiCult™ Mineral Oil

STERILE A

Document ID: FP09 I20 R01 B.3
Update: 08/05/2023

MATERIAL INCLUDED

Référence catalogue		
MINOIL050	1 x 50ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL100	1 x 100ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL500	1 x 500ml	FertiCult™ Mineral Oil (glass bottle)
MINOIL050P	1 x 50ml	FertiCult™ Mineral Oil (PETG bottle)
MINOIL100P	1 x 100ml	FertiCult™ Mineral Oil (PETG bottle)
MINOIL500P	1 x 500ml	FertiCult™ Mineral Oil (PETG bottle)

SUPPORT CLIENTS-SUPPORT TECHNIQUE

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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