FertiVit™ Cooling / FertiVit™ Warming

EN

FertiVit™ Cooling / FertiVit™ Warming

INTENDED USE

FertiVit™ Cooling and FertiVit™ Warming are a set of media for vitrification and warming of human oocytes and embryos (zygote till blastocyst).

For professional use only.

MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for 4-10 vitrification / warming procedures.

FertiVit™ Cooling kit (product code: FVC_KIT): contains one bottle each of the following media:

- Product code FPI005: 5 ml Pre-incubation medium ("PI")
- Product code FVC1001: 1 ml Cooling 1 ("C1")
- Product code FVC2001: 1 ml Cooling 2 ("C2")
- Product code FVC3001: 1 ml Cooling 3 ("C3")
- Product code FVC4001: 1 ml Cooling 4 ("C4")
- Product code FVC5001: 1.5ml Cooling 5 ("C5")

FertiVit™ Warming kit (product code: FWV_KIT): contains one bottle of each of the following media:

- Product code FWV1005: 7.5 ml Warming 1 ("W1")
- Product code FWV2001: 1 ml Warming 2 ("W2")
- Product code FWV3001: 1 ml Warming 3 ("W3")
- Product code FWV4001: 1 ml Warming 4 ("W4")
- Product code FWV5001: 1 ml Warming 5 ("W5")
- Product code FWV6001: 1 ml Warming 6 ("W6")

The media should be used in the order displayed above (the bottles may be in a different order in the box).

MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device, preferably closed system (FertiVit™ Cooling / FertiVit™ Warming was clinically validated with the use of Vitrisafe (Nextclincs) and High Security Vitrification device (HSV, CBS). Importantly, each laboratory should validate their own procedures when using FertiVit™ Cooling / FertiVit™ Warming in combination with other vitrification devices)

COMPOSITION

FertiVit™ Cooling and FertiVit™ Warming are ready to use. The media are HTF-based and contain HEPES and human serum albumin (12-20 g/liter depending on the individual medium; medicinal substance derived from human blood plasma). None of the media contain antibiotics.

FertiVit™ Cooling media have increasing concentrations of DMSO and Ethylene Glycol (EG). FertiVit™ Cooling 5 medium also contains ficol.

FertiVit™ Warming media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 – 7.50 (release criteria: 7.20-7.40)
- Osmolality (mOsm/kg):
 - Pre-incubation / Warming 6: 270-295 (release criteria: 270-290)
 - Warming 3: 805-865 (release criteria: 805-850)
 - Warming 4: 535-565
 - Warming 5: 405-435
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>; no growth
- Endotoxins test by Limulus Lysate methodology (USP <85>; < 0.25 EU/ml)
- Mouse embryo assay (one-cell mouse embryo development until blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 90%
- Use of Ph. Eur. or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
- seal of the container is opened or defect when the product is delivered.
- expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.

STORAGE/DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Stable after transport (max. 5 days) at elevated temperature (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification/warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm to room temperature (20-25 °C) before use. Warming 1 might be preheated to 37 °C before use.
- Open the necessary number of vitrification devices (for maximum volume load; check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of oocytes and embryos as described in [figure 1](#).

Notes:

- Do not use the same droplets/ media set-up for different patients.
- Set-up with 6-well dish (for cases with multiple embryos/oocytes): up to 5 vitrification cycles of the same patient

can be performed in the same media

- Set-up with droplets (for cases with few embryos/oocytes): the number of embryos / oocytes that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.
- Cooling 1 and 2 are only needed for vitrification of oocytes.
- Warming 5 is only needed for warming of zygotes and oocytes

Cooling protocol

- Expose oocytes / embryos sequentially to the following media as described in [figure 2](#).

- (*) Recommended pipetting technique in C5 ([figure 3](#)):
 - Place the oocyte / embryo on the bottom of the 1st drop of C5 (step A). Make sure the oocyte / embryo is in the drop. Preferably, do not extrude the complete volume of the C4. The oocyte / embryo will float.
 - Empty the pipette outside the drop (step B1).
 - Rinse the pipette with medium from drop 2 of C5 (step B2).
 - Aspirate a small volume of drop 2, aspirate the oocyte from drop 1, transfer to the bottom of drop 2 (step C).
 - Make sure the oocyte / embryo is in the drop. Flush in and out 2 times.
 - Place the oocyte / embryo in drop 3 and flush in and out 3 times until the shadow-like shape disappears (step D).
 - Now load the oocyte / embryo on the vitrification device.

Warming protocol

- Expose oocytes / embryos sequentially to the following media as described in [figure 4](#).

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to Fertipro and, if applicable, the Competent Authority of the EU Member State in which the user and/ or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for FertiVit™ Cooling and FertiVit™ Warming describes safety and performance characteristics for the media and is available on the website of Fertipro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact Fertipro NV for customer or technical support.

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