EN FR DE ES IT PT GR NL

GAIN™ medium

GAIN™ Medium

INDICATIONS FOR USE

GAIN™ medium is a single-step cell culture medium for use with human embryos and gametes. GAIN™ medium can additionally be used in all the procedures indicated in the figure below

- For semen washing and Intra Uterine Insemination (IUI).
- For oocytes handling/incubation in preparation of, or during fertilization by In Vitro Fertilization (IVF/ Intra Cytoplasmic Sperm Injection (ICSI)).
- For embryo culture from day 1 to expanded blastocyst stage.
- For embryo transfer (ET).

For professional use only.

GENERAL INFORMATION

The medium is complete and needs no further additives. GAIN™ medium has to be **pre-incubated in a CO₂ incubator with 5-6% CO₂ for at least 5 hours**, but **ideally, dishes should be prepared the day before use and incubated overnight** (with lid open) to obtain an optimal pH.

COMPOSITION

GAIN™ medium is a ready-to-use bicarbonate-buffered balanced salt solution, supplemented with 10 mg/l gentamicin (medicinal substance) and 3.5g/l human serum albumin (medicinal substance derived from human blood plasma) and phenol red.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20-7.45 (37 °C, 5% CO₂)
- Osmolality: 270-290 mOsm/kg
- Sterility test by the current Ph. Eur. 2.6.1/ USP <71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.25EU/ml
- One-cell mouse embryo assay (% blastocysts after 96 hours): ≥ 80%
- Human sperm survival assay (% motility compared with control, after 24 hours exposure to test medium at room temperature, no CO₂ incubation): ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com).

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes discoloured, cloudy, or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Keep in its original packaging until the day of use
- Product should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.

METHOD

Each laboratory should consult its own validated procedures, optimized for its individual medical program. When performing IUI or ET, follow the instructions of the specific catheters used.

Protocols described below are informative and provided as an example of good practice:

Fertilization and embryo culture (suggested protocol):

- Prepare dishes with six-eight 25-250µl droplets of GAIN™ medium covered with mineral oil (e.g. FertiCult™ Mineral Oil, FertiPro NV).
Note: for ICSI, it is convenient to prepare a dish with 2 central droplets (10µl droplet and a bigger droplet), and several peripheral droplets of GAIN™ medium.
- Before introducing oocytes / embryos in the droplet, incubate dishes overnight in a CO₂ incubator at 37 °C to ensure CO₂ and pH equilibration (see note below).
- After equilibration, place the oocytes in the droplets with GAIN™ medium.
A) For IVF: culture 1-4 trimmed OOCs (Cumulus Oocyte Complex) per droplet and inseminate with approximately 0.1 million sperm. Upon evaluation of pronuclei, carefully wash and transfer zygotes in fresh pre-equilibrated dishes with droplets of GAIN™ medium covered with mineral oil for further culture. (e.g. FertiCult™ Mineral Oil, FertiPro NV).
B) For ICSI: embryos can be placed per micro-droplet.
- For ICSI: Transfer oocytes to peripheral droplets of GAIN™ medium in ICSI dish (1 oocyte per droplet). Replace central 10µl of GAIN™ medium in ICSI dish by 10µl (prewarmed) PVP medium (e.g. 10% PVP in FertiCult™ Flushing medium, FertiPro NV). Transfer 1-2µl of prepared semen to central PVP drop in ICSI dish. Select spermatozoa in PVP drop for ICSI, wash selected sperm in the other central drop of GAIN™ medium, and subsequently perform ICSI of oocytes in the peripheral GAIN™ medium drops of the same dish. Further culture injected oocytes overnight in CO₂ incubator at 37 °C until embryo transfer.

- GAIN™ medium is a single-step medium designed for continuous embryo culture up to expanded blastocyst stage. However, if culture medium is refreshed, it is recommended to do this on day 2 or early day 3. In this case, make sure that the fresh medium is pre-incubated before transfer of the embryos.

- Transfer of embryos according to the instructions of the specific catheters used.

Notes:

Optimal pH and atmospheric pressure For optimal embryo culture, GAIN™ medium should be used at a pH of 7.28 ± 0.05. Different factors like the CO₂ concentration in the incubator and atmospheric pressure (which decreases at higher altitudes) have an effect on the pH after equilibration in the incubator¹. Therefore, we strongly advise to measure the pH under culture conditions with 5% CO₂ in the incubator in order to determine the CO₂ concentration setting of the incubator which result in an optimal pH of 7.28.

Washing of spermatozoa (suggested protocol)

Washing of spermatozoa can be done at room temperature or at 37 °C. (Pre-equilibration with CO₂ is not required)

- Add 5ml GAIN™ medium to the native semen sample and mix. Centrifuge for 5-10 minutes at approximately 300g.
- Remove supernatant and leave about 0.5ml of semen in the centrifuge tube.
- Add 5ml GAIN™ medium to the test tube. Mix the solution gently until the pellet is completely dissolved.
- Centrifuge again for 5-10 minutes at 300g

- Resuspend in a suitable volume of GAIN™ medium.

STORAGE/DISPOSAL INSTRUCTIONS

- Store GAIN™ medium between 2-8 °C.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from sunlight.
- The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic techniques should be used to avoid possible contamination, even though the products contain gentamicin.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/ removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for GAIN™ medium describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by downloading the following QR code:

For further questions regarding the to the safety and performance, please contact FertiPro NV. (see below: customer or technical support).

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