Hyaluronidase in FertiCult™ Flushing medium

Hyaluronidase in FertiCult™ Flushing medium

Hyaluronidase in FertiCult™ Flushing medium

INDICATIONS FOR USE

Hyaluronidase in FertiCult™ Flushing medium is used in the oocyte denudation process in preparation of intracytoplasmic sperm injection (ICSI) or other Assisted Reproductive Technologies. Hyaluronidase digests the hyaluronic acid between the cumulus cells, which makes it easier to remove the cumulus mechanically.

For professional use only.

GENERAL INFORMATION

The medium is complete and needs no further additives. The medium contains HEPES; **no CO2 incubation required**. Hyaluronidase in FertiCult™ Flushing medium does not contain any antibiotics

COMPOSITION

Hyaluronidase in FertiCult™ Flushing medium is a ready-to-use HEPES-buffered medium which also contains bicarbonate, physiologic salts, glucose, lactate, pyruvate, human serum albumin (4.0g/l, medicinal substance derived from human blood plasma) and 80 IU/ml hyaluronidase from bovine origin.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH criteria: 7.30-7.90 (release criteria: 7.30-7.60)
- Osmolality: 270-290 mOsm/kg
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.5EU/ml
- Sterility test by the current Ph. Eur. 2.6.1./ USP <71>: No growth
- One-cell mouse embryo assay (% blastocysts after 96 hours, exposure time to test medium: 30 minutes): ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
 - seal of the container is opened or defect when the product is delivered.
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- For product code HYA010: Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple opening/warming cycles of the medium. Discard excess (unused) media.
- Keep in its original packaging until the day of use

METHOD: example of suitable protocol

- Prepare a dish containing 1 droplet of Hyaluronidase in FertiCult™ Flushing medium (100µL) and 3-5 droplets (100µL) of an oocyte washing medium (e.g. FertiCult™ Flushing medium), all under mineral oil (e.g. FertiCult™ Mineral Oil) and warm to a temperature of 37 °C.
- Place COC's in the hyaluronidase in FertiCult™ Flushing medium (up to 5 oocytes maximum) for about 30 seconds.
- Using a fine glass pipette, transfer the partially denuded oocytes in the first washing droplet.
- Remove the corona cells by pipetting the oocytes.
- Use the other droplets to further wash the denuded oocytes.

Notes:

- Hyaluronidase in FertiCult™ Flushing medium is HEPES-buffered. Incubation in a CO2 incubator will lower the pH.
- Each laboratory should consult its own validated procedures, optimized for its individual media program.

STORAGE/DISPOSAL INSTRUCTIONS

- Store Hyaluronidase in FertiCult™ Flushing medium between 2-8 °C.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from sun/light.
- The products are stable after transport (max. 5 days) at elevated temperatures (≤ 37 °C).
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic techniques should be used to avoid possible contamination.
- Always wear protective clothing when handling specimens.
- All blood products should be treated as

potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non- elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.

- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- The hyaluronidase used in this product is derived from bovine testis and is certified with a Transmissible Spongiform Encephalopathy (TSE) risk evaluation Certificate of Suitability (CEP). The animals from which the hyaluronidase is derived, are determined "fit for human consumption" and originate from countries with "negligible BSE risk", as determined in Resolution "Recognition of the Bovine Spongiform Encephalopathy Risk Status of Member Countries", adopted by the World Organisation for Animal Health (OIE). According to the World Health Organization (WHO) guidelines on Tissue Infectivity Distribution in Transmissible Spongiform Encephalopathies (2010) and ISO22442-1, testes from bovine source are classified as "Tissues with no detected infectivity".
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV, and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for Hyaluronidase in FertiCult™ Flushing medium describes safety and performance characteristics for the media and is available on the website of FertiPro NV. (www.fertipro.com) or by downloading the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

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