Sil-Select Plus™

OUTDATED

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
LT	puslapis	11
NO	side	12
RO	pagină	13
SK	stran	14
TR	Sayfa	15
HU	oldal	16
DK	side	17

Sil-Select Plus™

STERILE A

Document ID: FP09 I13_Plus R01 E.9
Update: 28/03/2024

MATERIAL INCLUDED

Catalogue number		
SIP025LO	25x 2.5ml	Sil-Select Plus™ Lower layer
SIP025UP	25x 2.5ml	Sil-Select Plus™ Upper layer
SIP020LO	5x 20ml	Sil-Select Plus™ Lower layer
SIP020UP	5x 20ml	Sil-Select Plus™ Upper layer
SIP80_100	1x 100ml	Sil-Select Plus™ 80%

CUSTOMER-TECHNICAL SUPPORT

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

Procedure packs		
SIP008	8x 2.5ml	Sil-Select Plus™ Lower layer
	8x 2.5ml	Sil-Select Plus™ Upper layer
	8x 5ml	Sil-Select Plus™ Sperm Washing/Insemination medium
SIP016	2x 20ml	Sil-Select Plus™ Lower layer
	2x 20ml	Sil-Select Plus™ Upper layer
	5x 20ml	Sil-Select Plus™ Sperm Washing/Insemination medium
SIP050	1x 50ml	Sil-Select Plus™ Lower layer
	1x 50ml	Sil-Select Plus™ Upper layer
SIP100	1x 100ml	Sil-Select Plus™ Lower layer
	1x 100ml	Sil-Select Plus™ Upper layer

Single bottles of Sil-Select Plus and Sil-Select Plus Sperm Washing/Insemination medium are all CE marked (Notified Body number 2797) and can be ordered separately. Assembler of the procedure packs is FertiPro NV.

INTENDED USE AND COMPOSITION

Sil-Select Plus is a ready-to-use gradient system for semen preparation. The following gradients are available:

- Sil-Select Plus Upper layer (45%)
- Sil-Select Plus 80%
- Sil-Select Plus Lower layer (90%)

Sil-Select Plus gradients consist of silane-coated colloidal silica particles suspended in HEPES-buffered EBSS (Earle's balanced salt solution) supplemented with 0.4-2.2g/l human serum albumin (medicinal substance derived from human blood plasma). Some product codes are supplemented with 10mg/l gentamicin (medicinal substance). The Sil-Select Plus system can be used for semen preparation for Intra Uterine Insemination (IUI), In Vitro Fertilization (IVF) and IntraCytoplasmatic Sperm Injection (ICSI).

Sil-Select Plus Sperm Washing/Insemination medium is intended for use in assisted reproductive technologies, more specifically for sperm washing/holding and for the production of density gradients. The medium can be used for IUI. Sil-Select Plus Sperm Washing/Insemination medium is a HEPES-buffered medium which also contains bicarbonate, physiologic salts, glucose, lactate and human serum albumin (4.0g/l, medicinal substance derived from human blood plasma). The product is also available with phenol red and gentamicin (10mg/l, medicinal substance).

For professional use only.

PRODUCT SPECIFICATIONS

Sil-Select Plus Upper layer/Lower layer/80% media:

- Chemical composition
- pH criteria: 7.20-7.90 (release criteria: 7.20-7.60)
- Osmolality: 280-310 mOsm/kg (Upper layer)
290-330 mOsm/kg (Lower layer)
290-330 mOsm/kg (80%)
- Density: 1.0500-1.0700 (Upper layer)
1.1050-1.1150 g/ml (Lower layer)
1.097-1.107 g/ml (80%)
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.5EU/ml
- Sterility test according to the current Ph. Eur. 2.6.1/USP <71>: No growth
- Human sperm survival assay (% motility compared with control after 4 hours): ≥ 80%
- Human sperm survival assay (% motility compared with control after 24 hours): ≥ 75%
- Not MEA tested
- Use of Ph Eur or USP grade products if applicable
- A certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

Sil-Select Plus Sperm Washing/Insemination medium:

- Chemical composition
- pH criteria: 7.30-7.90 (release criteria: 7.30-7.60)
- Osmolality: 270-290 mOsm/kg
- Sterility test according to the current Ph. Eur. 2.6.1/USP<71>: No growth
- Endotoxin test by Limulus Amebocyte Lysate methodology (USP <85>): < 0.25EU/ml
- One-cell Mouse Embryo Assay (blastocysts after 96h): ≥ 80% after 60 min of exposure
- Human Sperm Survival Assay (% motility compared with control after 24 hours): ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- Certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

CALCULATIONS OF G-FORCES

The g-force of your centrifuge can be calculated using this formula:

or

$$g = 1.118 \times r \times \text{rpm}^2$$
$$\text{rpm} = \text{Square root } \{g / (1.118 \times r)\}$$

r = radius of centrifuge in mm
rpm = rotations per minute / 1000

Example 1
r = 100 mm
rpm = 1800 rotations per minute
g = 1.118 x 100 x 3.24 = 362 g

Example 2
r = 100 mm
g = 350g
rpm = SQR (350 / (1.118 x 100)) = 1.77
rpm = 1770 rotations per minute

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes discoloured (if medium contains phenol red), cloudy, or shows any evidence of microbial contamination;
 - seal of the container is opened or defect when the product is delivered;
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Product which include gentamicin should not be used on a patient that has a known allergy to gentamicin or similar antibiotics.
- Depending on the number of procedures that will be performed on one day, remove the required volume of medium under aseptic conditions in an appropriate sterile recipient. This is in order to avoid multiple openings/warming cycles of the medium. Discard excess (unused) media.
- Keep in its original packaging until the day of use.

METHOD

Importantly: each laboratory should consult its own optimized and validated laboratory procedures.

FRESH SEMEN SAMPLES

Mix the density gradient bottles by 5 bottle inversions before use.

- 1 Bring all components of the system and samples to room temperature or to 37 °C.
- 2 Transfer 2.5ml of Sil-Select Plus Upper layer into a sterile disposable centrifuge tube.
- 3 Place 2.5ml of Sil-Select Plus Lower layer under the Upper layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower layer. This two layer gradient is stable for up to two hours.
- 4 Gently place up to 2,5ml of liquefied semen onto the Upper layer using a transfer pipette. Do not use a higher volume than the volume of the individual gradient layers or more than 10⁹ cells.
- 5 Centrifuge for 15 to 18 minutes at 350g to 400g. When this centrifugation is completed you may not be able to visibly see a pellet. If so, it is essential to continue the procedure with a second centrifugation of 3 to 5 minutes.
- 6 Remove supernatant down to the pellet.
- 7 Add 2-5ml of sperm washing medium (e.g. FertiCult Flushing medium or Sil-Select Plus Sperm Washing/Insemination medium, FertiPro NV) and resuspend the pellet.
- 8 Centrifuge for 8 to 10 minutes at 300g. Higher sperm concentration will require the maximum 10 minutes centrifugation to ensure a complete and thorough sperm wash.
- 9 Remove supernatant down to the pellet and repeat steps 7 and 8.
- 10 Remove supernatant and replace with a suitable volume of appropriate medium.

FROZEN SEMEN SAMPLES

Mix the density gradient bottles by 5 bottle inversions before use.

- 1 Bring all components of the system and samples to room temperature or to 37 °C.
- 2 Transfer 1ml of Sil-Select Plus Upper layer into a sterile disposable centrifuge tube.
- 3 Place 1ml of Sil-Select Plus Lower layer under the Upper layer. Take care that the two layers are distinctly separated. This is done by placing the tip of the pipette on the bottom of the test tube and slowly dispensing the Sil-Select Plus Lower layer. This two-layer gradient is stable for up to two hours.

- 4 Gently place the thawed semen sample onto the Upper layer using a transfer pipette (0.5 - 1ml).
- 5 Centrifuge for 15-20 minutes at 350g.
- 6 Remove supernatant down to no less than the 0.5ml mark above the pellet.
- 7 Add 2-5ml of sperm washing medium (e.g. FertiCult Flushing medium or Sil-Select Plus Sperm Washing/Insemination medium, FertiPro NV) and resuspend the pellet.
- 8 Centrifuge for 8 to 10 minutes at 300g.
- 9 Remove supernatant down to the pellet and repeat steps 7 and 8.
- 10 Remove supernatant and replace with a suitable volume of appropriate medium.

STORAGE/DISPOSAL INSTRUCTIONS

- Store products without gentamicin between 2-25 °C. Once opened: store between 2-8 °C.
- Store products with gentamicin between 2-8 °C.
- The products can be used up to 7 days after opening when sterile conditions are maintained and the products are stored at 2-8 °C.
- Keep away from (sun)light
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

PRECAUTIONS

- Aseptic technique should be used to avoid possible contamination, even when the products contain gentamicin.
- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administrated, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven viruses transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro NV and, if applicable, to the competent authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for Sil-Select Plus and Sil-Select Plus Sperm Washing/Insemination medium describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by scanning the following QR codes:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1	
	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture
	Sterilized using aseptic processing techniques
	Temperature limit
	Single sterile barrier system
	Single sterile barrier system with protective packaging outside.
	Do not use if package is damaged and consult instructions for use.
	Do not resterilize
	Contains a medicinal substance
	Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745	
	CE marking by Notified Body 2797
Symbol as defined in 21 CFR 801.15	
	US Federal law restricts this device to sale by or on the order of a physician.
In compliance with MDR-article 22 on procedure packs, these CE marking is not applicable for procedure packs with product codes SIP008/SIP008PRG/SIP016/SIP016PRG/SIP050/SIP100. Procedure packs with different medical devices packed together cannot bear an additional CE marking themselves.	
Note	

Sil-Select Plus™

STERILE A

Document ID: FP09 I13_Plus R01 E.9
Update: 28/03/2024

MATÉRIEL INCLUS

Référence catalogue		
SIP025LO	25x 2.5ml	Sil-Select Plus™ Lower layer
SIP025UP	25x 2.5ml	Sil-Select Plus™ Upper layer
SIP020LO	5x 20ml	Sil-Select Plus™ Lower layer
SIP020UP	5x 20ml	Sil-Select Plus™ Upper layer
SIP80_100	1x 100ml	Sil-Select Plus™ 80%

SUPPORT CLIENTS-SUPPORT TECHNIQUE

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

Paquets de procédures		
SIP008	8x 2.5ml 8x 2.5ml 8x 5ml	Sil-Select Plus™ Lower layer Sil-Select Plus™ Upper layer Sil-Select Plus™ Sperm Washing/ Insemination medium
SIP016	2x 20ml 2x 20ml 5x 20ml	Sil-Select Plus™ Lower layer Sil-Select Plus™ Upper layer Sil-Select Plus™ Sperm Washing/ Insemination medium
SIP050	1x 50ml 1x 50ml	Sil-Select Plus™ Lower layer Sil-Select Plus™ Upper layer
SIP100	1x 100ml 1x 100ml	Sil-Select Plus™ Lower layer Sil-Select Plus™ Upper layer

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