VitriFreeze ES™

VitriThaw ES™

OUTDATED

Available languages

EN	page	2
FR	page	3
DE	Seite	4
ES	página	5
IT	pagina	6
PT	página	7
GR	σελίδα	8
NL	pagina	9
BG	страница	10
RO	pagină	11
TR	Sayfa	12
LT	puslapis	13
NO	side	14
DK	side	15
HU	oldal	16
CZ	stránka	17

VitriFreeze ES™

VitriThaw ES™

Document ID: FP09 I46 02 R01 E.3
Update: 15/04/2024

MATERIAL INCLUDED

Catalogue number

VF_KIT1_ES	VitriFreeze ES - kit
VT_KIT1_ES	VitriThaw ES - kit

CUSTOMER-TECHNICAL SUPPORT

FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

INTENDED USE

VitriFreeze ES™ and VitriThaw ES™ are a set of media for vitrification and thawing of human embryos (zygote till expanded blastocyst stage).

For professional use only.

MATERIAL INCLUDED WITH THE KIT

One kit provides sufficient medium for a maximum of 6 vitrification/ warming procedures, depending on the protocol used.

VitriFreeze ES™ kit (product code VF_KIT1_ES): contains one bottle of the following media:

- Product code VP1E005: 5 ml VitriFreeze ES Pre-incubation medium ("PI")
- Product code VFE1001: 1 ml VitriFreeze ES Freezing medium 1 ("F1")
- Product code VFE2001: 1 ml VitriFreeze ES Freezing medium 2 ("F2")
- Product code VFE3001: 1 ml of VitriFreeze ES Freezing medium 3 ("F3")

VitriThaw ES™ kit (product code: VT_KIT1_ES): contains one bottle of the following media:

- Product code VTE1005: 5 ml of VitriThaw ES Thawing medium 1 ("T1")
- Product code VTE2001: 3.2 ml of VitriThaw ES Thawing medium 2 ("T2")
- Product code VTE3001: 1 ml of VitriThaw ES Thawing medium 3 ("T3")
- Product code VTE4001: 1 ml of VitriThaw ES Thawing medium 4 ("T4")
- Product code VTE5001: 1 ml of VitriThaw ES Thawing medium 5 ("T5")

The media should be used in the order displayed above (the bottles may be in a different order in the kit).

MATERIAL NOT INCLUDED BUT NEEDED

- Freezing tank with liquid nitrogen
- Vitrification device, preferably closed system (VitriFreeze ES™/ VitriThaw ES™ were clinically validated with use of the following closed devices: High Security Vitrification straw (HSV; CBS) or Vitrisafe (Nextclinics). Importantly, each laboratory should validate their own procedures when using VitriFreeze ES™/VitriThaw ES™ in combination with other vitrification devices)

COMPOSITION

VitriFreeze ES™ and VitriThaw ES™ are ready to use. The media contain PBS and human serum albumin (11-20 g/liter depending on the individual medium, medicinal substance derived from human blood plasma). None of the media contain antibiotics.

VitriFreeze ES™ media have increasing concentrations of DMSO and Ethylene glycol (EG). VitriFreeze ES™ Freezing 3 also contains sucrose and ficoll. VitriThaw ES™ media have decreasing concentrations of sucrose.

PRODUCT SPECIFICATIONS

- Chemical composition
- pH: 7.20 – 7.40
- Osmolality (mOsm/kg):
 - VitriFreeze ES™Pre-incubation medium: 270-295 (release criteria: 270-290)
 - VitriThaw ES™Thawing medium 2: 805-865 (release criteria: 805-850)
 - VitriThaw ES™Thawing medium 3: 535-565
 - VitriThaw ES™Thawing medium 4: 405-435
- Sterility test by the current Ph. Eur. 2.6.1. / USP <71>: no growth
- Endotoxins test by Limulus Lysate (LAL) methodology (USP <85>): < 0.25 EU/ml

- Mouse Embryo Assay (one-cell mouse embryo development until blastocysts stage after 96 hours following a sequential exposure to the cooling and warming media) ≥ 80%
- Use of Ph Eur or USP grade products if applicable
- The certificate of analysis and MSDS are available upon request or can be downloaded from our website (www.fertipro.com)

WARNINGS BEFORE USE

- Do not use the product if:
 - it becomes cloudy, or shows any evidence of microbial contamination.
 - seal of the container is opened or defect when the product is delivered.
 - expiry date has been exceeded.
- Do not freeze before use.
- Do not re-sterilize after opening.
- Freezing medium 2 may contain small salt precipitates which do not have an impact on product performance/safety.

STORAGE /DISPOSAL INSTRUCTIONS

- Store between 2-8 °C.
- Keep away from (sun)light.
- The products can be used up to 7 days after opening, when sterile conditions are maintained and the products are stored at 2-8 °C.
- The devices need to be disposed in accordance with local regulations for disposal of medical devices.

METHOD

We strongly advice to read through all the steps of the vitrification/ warming procedure before starting the procedure.

Preliminary steps

- Mix all bottles of the kit well and warm the media to room temperature (20-25 °C).
- Open the necessary number of vitrification devices, taking into account that 1 device can hold 1-2 embryos, (for maximum volume load: check the instructions for use of the device you are using). Conveniently place the separate parts of the vitrification device on the workbench for easy access later in the procedure.
- Prepare the following dishes for vitrification or warming of embryos as described in **figure 1**.

Notes:

- Do not use the same droplets/ media set-up for different patients.
- Set-up with 4-well dish (for cases with multiple embryos): up to 5 vitrification cycles of the same patient can be performed in the same media.
- Set-up with droplets (for cases with few embryos): the number of embryos that can be processed in the same medium is limited. Droplets without oil overlay should be prepared immediately before processing.

Vitrification protocol

- Expose the embryos sequentially to the following media as described in **figure 2**.

- (*) Recommended pipetting technique for step using F3 (**figure 3**):
 - Place the embryo on the bottom of the 1st drop of F3 (step A); Make sure the embryo is in the drop. Preferably, do not extrude the complete volume of the F2. The embryo will float.
 - Empty the pipette outside the drop (step B1).
 - Rinse the pipette with medium from drop 2 of F3 (step B2).
 - Aspirate a small volume of drop 2, aspirate the embryo from drop 1, transfer to the bottom of drop 2 (step C).
 - Make sure the embryo is in the drop. Flush in and out 2 times.
 - Place the embryo in drop 3 and flush in and out 3 times until the shadow-like shape disappears (step D).
 - Now load the embryo on the vitrification device.

Thawing protocol

- Expose the embryos sequentially to the following media as defined in **figure 4**.

- To allow a more transient removal of cryoprotectants (highly recommended for early embryo stages), an additional thawing step can be added between the T1 and T2 step, in which embryos are exposed for 1 minute to a 1:1 mixture T1/T2.

PRECAUTIONS

- All blood products should be treated as potentially infectious. Source material used to manufacture this product was tested and found non-reactive for HbsAg and negative for Anti-HIV-1/-2, HIV-1, HBV, and HCV. Furthermore, source material has been tested for parvovirus B19 and found to be non-elevated. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents.
- Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.
- Always wear protective clothing when handling specimens.
- Aseptic techniques should be used to avoid possible contamination.
- Any serious incident (as defined in European Medical Device Regulation 2017/745) that has occurred should be reported to FertiPro and, if applicable, the Competent Authority of the EU Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP for VitriFreeze ES™ and VitriThaw ES™ describes safety and performance characteristics for the media and is available on the website of FertiPro NV (www.fertipro.com) or by using the following QR code:

For further questions regarding to the safety and performance, please contact FertiPro NV for customer or technical support.

SYMBOLS GLOSSARY

Symbols as defined in ISO 15223-1

REF	Catalogue number
LOT	Batch code
	Use-by date
	Manufacturer
	Keep away from sunlight
MD	Medical device
	Consult instructions for use
	Cautions
	Date of manufacture

STERILE A Sterilized using aseptic processing techniques

°C Temperature limit

Single sterile barrier system

Single sterile barrier system with protective packaging outside.

Do not use if package is damaged and consult instructions for use.

Do not resterilize

Contains a medicinal substance

Contains human blood or plasma derivatives

Symbol as defined in EU Commission MDR 2017/745

CE marking by Notified Body 2797

Symbol as defined in 21 CFR 801.15

US Federal law restricts this device to sale by or on the order of a physician.


Figure 2

stage	PI	F1	F2	F3* (see figure 3)	
zygote till expanded blastocyst	2 min	5-7 min	4-5 min	40-60 sec	*60 seconds time to: 1. place the embryo in F3 2. load on the vitrification device 3. insert in the outer straw 4. seal 5. plunge into the liquid nitrogen
DMSO / EG (%)	-	5	10	20	


Figure 4

	T1 (can be at 37 °C) gently swirl the straw in the medium to ensure more homogenous temperature	T2	T3	T4	T5	Pre-equilibrated culture medium	
zygote till expanded blastocyst	1 min	1-2 min	2-4 min	2-4 min	1-2 min	wash	cultivate
sucrose (M)	1	0.5	0.25	0.125	0		

VitriFreeze ES™

VitriThaw ES™

Document ID: FP09 I46 02 R01 E.3
Update: 15/04/2024

MATÉRIEL INCLUS

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VF_KIT1_ES	VitriFreeze ES - kit
VT_KIT1_ES	VitriThaw ES - kit

SUPPORT CLIENTS-SUPPORT TECHNIQUE

 FertiPro NV
Industriepark Noord 32
8730 Beernem / Belgium
Tel +32 (0)50 79 18 05
Fax +32 (0)50 79 17 99
URL: www.fertipro.com
E-mail: info@fertipro.com

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