

To whom it concerns,

FertiPro N.V. updated the instructions for use of SpermFreeze. A summary of the changes is listed in the table below.

Instructions for use	SpermFreeze
Previous version	FP09 I11 R01 D.1
New version	FP09 I11 R01 D.2
Changes	<ul style="list-style-type: none">- Some small editorial changes (symbol sperm handling is added, the bibliography and contact address are added to the black section, the wording "sterile filtration" has been replaced by "aseptic filtration"; the HSA concentration has not changed in the product, but is rounded in the IFU to 4 g/L, "only for the intended use" is deleted)- Wash medium is added in the section "Material not included"- Small correction in the product specifications: the sperm survival test is performed on density selected spermatozoa.- The temperature of the tap water used in step 2 of the thawing protocol is added: "<i>Place the straws in tap water for 5 minutes (room temperature or 37°C).</i>"- The Portuguese, Dutch and Greek translation are added to the instructions for use.

Name: Ellen Van Vlierberghe

Function: Regulatory Affairs Associate

Signature and stamp:


Ellen Van Vlierberghe
Regulatory Affairs Associate
FertiPro N.V.