


To whom it concerns,

FertiPro NV updated the instructions for use of Sil-Select Stock. A summary of the changes is listed in the table below.

Instructions for use	Sil-Select Stock
Previous version	FP09 I13_Stock R01 E.4
New version	FP09 I13_Stock R01 E.5
Changes	<ul style="list-style-type: none"> - Correction of available product codes in 'Material included': - Indications for use, intended user and target patient population are more clearly defined as separate items (regulatory requirement). - Small update in section 'Method': <ul style="list-style-type: none"> ➤ Use of syringe was removed ➤ Add 2-5ml of sperm washing medium to resuspend the pellet - Section 'Storage/disposal instructions': Store products between 2-8°C. - Statement on stability after transport removed - Update in symbol glossary - Addition of QR code: 'Other languages can be downloaded on our website' - General corrections in translations

Name: Aily Standaert
Function: Regulatory Affairs Associate
Date: 27/05/2025
Signature and stamp:


Aily Standaert
Regulatory Affairs
FertiPro N.V.

COMPANY ADDRESS

FertiPro nv
Industriepark Noord 32
B-8730 Beernem / Belgium

tel. +32 (0)50 79 18 05
fax +32 (0)50 79 17 99
info@fertipro.com
www.fertipro.com

RPR Brugge
VAT BE 0448 767 035
IBAN BE 98 7374 4024 9293
BIC KREDBEBB