

06/06/2025

To whom it may concern,

FertiPro NV updated the instructions for use of SpermMar Test IgA. A summary of the changes is listed in the table below.


Instruction for use	SpermMar Test IgA
Previous version	FP09 I01 R01 F.2
New version	FP09 I01 R01 F.3
Changes	- Make differentiation between Single and Complete kit more clear - Vortex latex particles after adjustment to room temperature - Removal of '5 times' in the protocol. - Update of the QR code to include the new instruction video on SpermMar Test IgA. - All other languages than English: general update in layout.

Name: Eline De Reu

Function: Regulatory Affairs Associate

Date: 06/06/2025

Signature and stamp:


Eline De Reu
Regulatory Affairs
FertiPro N.V.

COMPANY ADDRESS

FertiPro nv
Industriepark Noord 32
B-8730 Beernem / Belgium

tel. +32 (0)50 79 18 05
fax +32 (0)50 79 17 99
info@fertipro.com
www.fertipro.com

RPR Brugge
VAT BE 0448 767 035
IBAN BE 98 7374 4024 9293
BIC KREDBEBB