

LeucoScreen Plus

Improved leucocyte staining method

PRODUCT

LeucoScreen Plus is an improved histochemical staining kit for the determination of peroxidase-positive white blood cells in semen. The kit provides useful information in semen samples where the concentration of round cells exceeds 1 million/ml. LeucoScreen Plus has an equal capability to distinguish between peroxidase-positive and -negative round cells as the LeucoScreen kit (same assay performance), but has the advantage that the substrate is not carcinogenic. The color of the substrate in peroxidase-positive white blood cells is bluish-grey to black instead of brown (LeucoScreen). The kit has been designed for 40 days of analysis.

PRODUCT INSTRUCTION VIDEO

An instruction film on the use of LeucoScreen Plus is available on YouTube.



STORAGE

+25°C
+2°C
Store between 2-25°C

12 months
from date of manufacture

REGULATORY

Europe: CE-marked – Other regions: information available upon request

PRODUCT ORDER CODE

LEUCO_PLUS LeucoScreen Plus – 40 days of analysis